MICRO-MOTORS

User's manual • Manuel d'utilisation

Manual de uso • Benutzerhandbuch

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ENGLISH

TYPE

IS, AV, IML electric dental micromotors with brushes. Non-sterilizable, protected to the oil of the hand-pieces. Removable, sterilizable sleeve.

Interchangeable brushes.

DESCRIPTION

IS: without light, internal spray.

AV: micromotor with higher torque, self-ventilating.

IML: micromotor with internal spray and light with bulb and switch integral with the motor.

INTENDED USE

Product intended for professional use only. Use in dentistry for prophylaxis, general dentistry and endodontic work. Any use other than that for which this device is intended is prohibited and may prove dangerous.

TECHNICAL DATA

Classification: Class IIa in accordance with European Directive 93/42/EEC concerning

medical devices. This medical device is in compliance with the legislation in force.

Coupling: The most commonly used coupling in the world as per ISO Standard 3964.

Dimensions : AV: 24 x 92.6mm - IML: 23.3 x 86mm - IS: 23.5 x 81mm

Weight: AV: 111g - IML: 97g - IS: 75g

Speed: Maximal 40.000rpm.

Recommended rotation speed: 60rpm (minimum is variable, depending on the electronic power supply used) up 40.000rpm \pm 10%; rotation in both directions.

Nominal voltage: In accordance with ISO 11498 : 0-24Vdc

Continuous operation: AV: 1.5A; approx. 0.8N.cm - IML & IS: 1A; approx. 0.6N.cm

Intermittent operation: According to IEC 34-1, Type S3. 3A for 60s; 5A for 10s; air cooling in 3min or without air in 20min

Mechanical power: AV: 74.2W - IML: 54W - IS: 56W

Torque : AV: 2.1N.cm - IML & IS: 2.3N.cm

Idle current: 0.3A - 24Vdc

Current at max. power: AV: 5A - IML & IS: 4A

Internal resistance: AV: 2.15 Ω - IML: 2.61 Ω - IS: 2.73 Ω

Noise: In accordance with ISO 7785-1, less

than 62dBA at 45cm (17.72in)

Bulb: IML: 3.5V; 0.74A

Comments: The performance of the micromotor varies, depending on the type of electronic power supply used. We recommend that you use Satelec electronic control units for optimum micromotor performance.

Electromagnetic compatibility: Corresponds to the electromagnetic compatibility in accordance with IEC 60601-1-2. Refer to the instructions for the device with which the motor is used.

Electrical safety: According to IEC 60601-1 Standard (General safety for Medical Electrical Equipment). Refer to the instructions for the device with which the motor is used.

Information: The technical specifications, illustrations and dimensions contained in these instructions are given only as a guide. They may not be the subject of any claim. The manufacturer reserves the right to make technical improvements to its equipment, without amending these instructions. For all additional information, please contact Satelec at the address indicated on the back cover.

ASSEMBLY

Connecting the hoses: Check the rear of the motor and the joint on the hose are clean. Line up the connection components on the motor with those of the plug and insert until stop position.

Fully tighten tube retainer. Replace sleeve fig. 2-4.

Cooling: Via compressed air from the unit. Place the flow-meter on the coupling and regulate the flow at 8 normliter/min. for the IML micromotor, and at 3 normliter/min. for the IS micromotor fig. 1, except for the AV micromotor, which is self-ventilating.

Important: Liquid or spray must not be allowed to penetrate inside the micromotor, due to risk of irreversible damage.

Never connect an instrument with a micromotor that is running.

Changing the bulb: IML fig. 5. Box of 5 bulbs. Wear rubber gloves when carrying out this changing operation.

Changing the brushes: **Note:** fine work; keep the brushes dry, NEVER in contact with the lubricant, otherwise there is a risk of irreversible damage. Before removing the two carbon brushes, wipe the back of the

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To change the two carbon brushes, remove the used carbon brushes. Blow dry, clean air across the apertures at the rear of the body of the micromotor fig. 9. To reassemble, put the brushes in place and secure them fig. 10 with the spring.

Replace the sleeve: Take care not to damage the O ring when replacing the sleeve. IS fig. 7 - AV fig. 8 - IML fig. 6

Seals replacement: To change the O rings fig. 11. Only valid for the IS micromotor.

Changing of the O ring fig. 12, not valid for the IML micromotor.

Changing of the O ring fig 8 + 13, not valid for the AV micromotors.

MAINTENANCE

Important: Never lubricate the electric brush motor.

Only use maintenance components from Satelec. The use of other components can void the guarantee.

Precautions for maintenance: The universal precautions, in particular wearing of personal

protective equipment (gloves, goggles, etc.), should be complied with by medical personnel using or performing maintenance on medical devices that contaminated or potentially contaminated. Pointed and sharp instruments should be handled with great caution.

Suitable detergent: For manual cleaning-disinfection:

- Detergent or detergent-disinfectant (pH 6-9.5) recommended for cleaning-disinfection of dental or surgical instruments. Quaternary ammonium- and/or enzymebased surfactants.
- Do not use detergents that are corrosive or contain chlorine, acetone aldehydes or bleaches.

Manual cleaning-disinfection only: Attention: Electric motors with brushes are not suitable for automatic cleaning-disinfection (in a washer-disinfector machine).

Clean and disinfect using a clean cloth soaked in a suitable product.

- Never immerse in disinfectant solutions.
- Do not soak in physiological liquid (NaCl).
- Do not immerse in an ultrasonic bath.

Sleeve sterilisation: Only for IML micromotor (fig. 6).

Check that the steriliser and the water that is used are clean. After each sterilisation cycle, remove the device from the sterilising apparatus immediately, in order to reduce the risk of corrosion.

Conditioning process: Wrap the removal sleeve in approved packaging for steam sterilisation.

Important: The quality of sterilisation depend very much on the cleanliness of the device. Only perfectly clean devices may be sterilised. Do not use a sterilisation procedure other than the one described below:

- With steam after fractioned initial vacuum phase acc. to EN 285 or with Class B cycle acc. to EN 13060. Nominal temperature: 134°C (273°F).
- Duration: 3 or 18 min. both procedures have been validated according to ISO 17664.

Overhaul: Never disassemble the device. For any modification and repair, we recommend that you contact your regular supplier or Satelec directly. Satelec asks the user to have its dynamic instruments checked or inspected at least once a year.

Transport and storage conditions:
Temperature between -40°C (-40°F) and 70°C (158°F), relative humidity between 10% and 100%, atmospheric pressure 500hPa to 1060 hPa (7.3 to 15.3psi).

OTHER PRECAUTIONS FOR USE

The device must be used by a qualified person in accordance with the current legal provisions concerning industrial safety, health and accident prevention measures, and these working instructions. In accordance with these requirements, the operators:

- must only use operating devices that are in perfect working order; in the event of irregular functioning, excessive vibration, abnormal heating or other signs indicating malfunction of the device, the work must be stopped immediately; in this case, contact a repair centre that is approved by Satelec;
- must ensure that the device is used only for the purpose for which it is intended, must protect themselves, their patients and third parties from any danger, and must avoid contamination through the use of the product.

Rest the device on a suitable support to avoid risks of infection for yourself, the patient or third parties.

The device is not authorised for use in an explosive atmosphere (anaesthetic gas).

It is essential to use dry, purified compressed air in order to ensure the long working life of the device. Maintain the quality of the air and the water by regular maintenance of the compressor and the filtration systems. The use of unfiltered hard water will lead to early blockage of the tubes, connectors and spray cones.

DISPOSAL

 This device must be recycled.
 Electrical and electronic equipment may contain dangerous substances which constitute health and environmental hazards. The user must return the device to its dealer or establish direct contact with an approved body for treatment and recovery of this type of equipment (European Directive 2002/96/ EC).

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