

# DENTALAIRE

Meeting and Exceeding the Needs of Veterinary Dentistry since 1988  
www.DentalaireProducts.com  
800-866-6881

## February is Dental Month, Do You Need to Replace Your Dentalaire Scaling Tip?

*Air Sourced Dental Stations*

*Portable Dental Stations*

*SilentAire Compressors*

*SilentAire Oil-Less Compressors*

*Hand Instruments*

*X-Ray Generators*

*Digital Dental X-Ray Systems*

*Natural Bone Graft Material*

*Restorative Products*

*Endodontic Equipment | Supplies*


*Consumable Supplies*

*Much More*



### Check your Piezo Tips!

Simply hold your tip up to the corresponding illustration – if yours is worn down more than 2mm it is time to replace it!



Purchase a #39 Universal Tip and a #37 Interproximal tip and receive a #31 tip **FREE!**

- NEW!** - The #31 tip
- This scaling tip is used for lingual, buccal and approximal supragingival scaling
  - Recommended for removal of gross calculus

To receive your free #31 tip, please fax your Scaler Tip invoice (#39 and #37) to 714-540-9947. Offer ends March 31st.

**Did you know that every mm  
of tip wear reduces efficiency by 25%?**



Everything You Need in Veterinary Dentistry . . . . .