

# EQUIPMENT COMPARISON CHART

| Product # | Item                        | Hi Speed | Low Speed | Water & Air | Piezo Scaler | Air Scaler | Mobile Cart | Counter Top | Wall Mount | Table Mount | Fiber-optics | P.Button Turbine | Compressor |
|-----------|-----------------------------|----------|-----------|-------------|--------------|------------|-------------|-------------|------------|-------------|--------------|------------------|------------|
| DTP00494  | Dentalaire II               | X        | X         | X           |              |            | X           |             |            |             |              |                  |            |
| DTP00495  | Dentalaire II               | X        | X         | X           |              |            | X           |             |            |             |              |                  | X          |
| DTP00500  | Dentalaire II               | X        | X         | X           |              |            |             | X           |            |             |              |                  | X          |
| DTP00504  | Dentalaire III              | X        | X         | X           | X            |            |             | X           |            |             |              | X                | X          |
| DTP00503  | Dentalaire III with Cabinet | X        | X         | X           | X            |            |             | X           |            |             |              | X                | X          |
| DTP00505  | Dentalaire III              | X        | X         | X           |              | X          |             | X           |            |             |              |                  | X          |
| DTP00506  | Dentalaire III              | X        | X         | X           |              | X          |             | X           |            |             | X            | X                | X          |
| DTP00507  | Dentalaire III              | X        | X         | X           |              | X          | X           |             |            |             |              |                  | X          |
| DTP00508  | Dentalaire III              | X        | X         | X           |              | X          | X           |             |            |             |              |                  |            |
| DTP00509  | Ultimate Vet Dental Station | X        | X         | X           | X            |            | X           |             |            |             | X            | X                | X          |
| DTP00510  | Dentalaire III              | X        | X         | X           | X            |            | X           |             |            |             | X            | X                |            |
| DTP00550  | Cabinet to Countertop Unit  |          |           |             |              |            |             |             |            |             |              |                  |            |
| DTP00600  | Two Handpiece               | X        | X         | X           |              |            |             |             | X          |             |              |                  |            |
| DTP00601  | Two Handpiece               | X        | X         | X           |              |            |             |             | X          |             |              |                  | X          |
| DTP00605  | Three Handpiece             | X        | X         | X           |              | X          |             |             | X          |             |              |                  |            |
| DTP00606  | Three Handpiece             | X        | X         | X           |              | X          |             |             | X          |             |              |                  | X          |
| DTP00607  | Dentalaire Basic            | X        | X         | X           |              |            |             |             |            | X           |              |                  |            |
| DTP00608  | Dentalaire Basic            | X        | X         | X           |              |            |             |             |            | X           |              |                  | X          |

## SCALERS / MICROMOTORS

|          |  |
|----------|--|
| DTP00772 | Scalex "Classic" Ultrasonic Scaler: Has Cavitron Inserts, 25K        |
| DTP00787 | Scalex "Classic" Ultrasonic Scaler: Has Cavitron Inserts, 30K        |
| DTP103RM | Ram 1000 Complete Prophy and Restorative Unit, 35,000 RPM Micromotor |
| DTP00780 | Rotex 780 Micro Motor, 30,000 RPM                                    |
| DTP00773 | Pressurized Water Bottle System                                      |
| DTP00776 | Scalex 880 Scaler, Piezo Electric Scaler                             |
| DTP00770 | Amdent US30 Scaler, Piezo Electric Scaler                            |
| DTP00789 | MaxiSonic Scaler Polisher Combo; Micromotor / Piezo Electric Scaler  |
| DTP00769 | MaxiSonic S, Piezo Electric Scaler                                   |

All **Countertop Dental Systems** come with a **1000ml. water bottle**, which mounts on the inside of the compressor.

The **Mobile, Wall Mount, and Table Mount Systems** come with a quick connect **1500ml. water bottle**.

**Distilled water is RECOMMENDED!**