


FRICITION GRIP BURS

Many types of burs are available. The cross-cut fissure bur (699, 700, 701, 701L, 702) is a good all-purpose bur that can be used to section teeth and gain endodontic access. It is also useful as a general purpose bur.


The round bur (1/2, 1, 2, 3, 4, 6, 8) is used to gain endodontic access, to prepare cavities, and to mark the cutting depth in crown preparation. The inverted-cone bur (33-1/2, 34, 35, 37) is used for undercutting cavity preparation to create a mechanical interlock. The pear-shaped bur (330, 331, 332) is used for root canal access and for undercutting dentin.

ROUND




- DTP10046 – Round Bur #1/2, Pk. of 5
- DTP10047 – Round Bur #1, Pk. of 5
- DTP10048 – Round Bur #2, Pk. of 5
- DTP10050 – Round Bur #3, Pk. of 5
- DTP10052 – Round Bur #4, Pk. of 5
- DTP10054 – Round Bur #6, Pk. of 5
- DTP10056 – Round Bur #8, Pk. of 5

INVERTED CONE




- DTP10033 – Cone Shape #33-1/2, Pk. of 5
- DTP10034 – Cone Shape #34, Pk. of 5
- DTP10035 – Cone Shape #35, Pk. of 5
- DTP10037 – Cone Shape #37, Pk. of 5

PEAR SHAPE



- DTP10017 – Pear Shape #330, Pk. of 5
- DTP10018 – Pear Shape #331, Pk. of 5
- DTP10020 – Pear Shape #332, Pk. of 5

CROSSCUT TAPER FISSURE



- DTP10025 – Cross Cut #699, Pk. of 5
- DTP10026 – Cross Cut #700, Pk. of 5
- DTP10028 – Cross Cut #701, Pk. of 5
- DTP10030 – Cross Cut #701L, Pk. of 5
- DTP10032 – Cross Cut #702, Pk. of 5
- DTP10021 – Cross Cut Fissure #557, Pk. of 5
- DTP10019 – Cross Cut Fissure #558, Pk. of 5



MAGNETIC BUR BLOCK ROUND

DTP10205

REPLACEMENT TIPS

FOR ULTRASONIC, MEGASONIC, and all AIR SCALERS

- DTP10880 – Cavitron Insert TFI-10 25K (Fine)
- DTP10881 – Cavitron Insert TFI-10 30K (Fine)
- DTP10879 – Cavitron Insert (Beavertail) TFI-9 25K (Wide)
- DTP100UN – Megasonic Tip (Universal)
- DTP100FL – Megasonic Tip (Flat)
- DTP100PR – Megasonic Tip (Perio)
- DTP100CT – CatheterTip f/ Megasonic
- DTP100TW – Tip Wrench f/ Megasonic
- DTP10876 – Star Titan Universal Tip
- DTP10877 – Star Titan Perio Tip
- DTP10878 – Star Titan Rebuild Kit
- DTP10882 – Sonix, Lynx, Airstream (Universal)
- DTP10883 – Sonix, Lynx, Airstream (Perio)
- DTP10884 – Sonix, Lynx, Airstream (Sickle)
- DTP10870 – Tip for Air Scaler (Sickle)
- DTP10875 – Tip for Air Scaler (Universal)
- DTP00576 – Replacement Spindle for Air Scaler
- DTP10872 – Tip for Air Scaler (Perio)
- DTP10885 – Sonus V Stack with Nosecone
- DTP10886 – Universal Tip to Sonus V Stack
- DTP10887 – Flat Tip to Sonus V Stack
- DTP00784 – Universal Tip #39
- DTP00785 – Flat Tip #32
- DTP00786 – Special / Perio Tip #33
- DTP00791 – Interproximal #37

TIP DIAGRAM

