

# ***The Ultimate Veterinary Dental Station***

*The only unit you will need for all of your dental procedures . . . With the DentaIaire "Ultimate" Dental Station there are no "options" needed. The unit is totally complete with everything you need to perform all of your dental procedures.*

## **Features**

Virtually soundless and powerful air compressor to run all day, everyday for maximum productivity

Two handpiece automatic control system

Tank powder coated inside to prevent corrosion

Dual filtration system, which gives you cleaner air

Three way (air, water and spray) syringe

Built-in water system with 1.5 liter water bottle

Built-in ultrasonic scaler with autoclavable tips

High speed fiberoptic handpiece with push button turbine

Low speed polishing handpiece with prophylaxis angle prophylaxis paste and cups

Telescopic height adjustment from 38" down to 24" for convenient storage


Work area 22" x 18" with stainless steel removable tray

Larger 9 outlet multiple strip for extra accessories



DTP00509

## **Special Features Include:**



### **High Speed Fiberoptic Handpiece**

with oral illumination system. Push button control for ease of bur changing.



### **Built-in Ultrasonic Piezo Electric Scaler**

with Universal and Flat Scaling Tips. Runs cool and quiet at 30,000 cycles per second.



### **Low Speed Handpiece**

with "reciprocating" prophylaxis angles. No more heat – no more splatter and no more tangles.

**For additional information on DentaIaire products contact your local distributor or call DENTALAIRE direct at (800) 866-6881**

# **DENTALAIRE™**

17150 Newhope Street  
Suite 407  
Fountain Valley, CA 92708  
E-MAIL: DENTALAIRE@aol.com