



DENTALAIRE

Meeting and exceeding the needs of veterinary dentistry since 1988.

2026

DIGITAL CATALOG



VETERINARY DENTAL EQUIPMENT



Meeting and exceeding the needs of veterinary dentistry since 1988.

At Dentalaire, we understand the importance of being able to provide quality veterinary care, especially when it comes to an animal's dental health.

With that in mind, it is **our vision to be the preeminent dental supplier in the veterinary market**. By offering new and innovative products and partnering with key dental professionals across the country, we are quickly nearing this goal.

From dental stations to imaging systems, we strive to provide veterinarians and dental professionals with the resources and equipment needed to help grow and improve their veterinary dentistry department, while also exceeding industry standards.

Visit our website to learn more: www.dentalaire.com

Dentalaire Support

Business Hours:

Monday - Friday: 7:30am - 4pm (PST)

Saturday - Sunday: Closed

Dentalaire Customer Service and Tech Support

Toll Free: (800)866-6881

E-mail: customerservice@dentalaire.com

Address:

17742 Mitchell North, STE B,
Irvine, CA, 92614

Website:

www.dentalaire.com



Scan the QR code to explore the most up to date Dentalaire catalog on your mobile device.

Table of Contents

Last Revision 2.11.2026

DTX Imaging System	6
DTX Imaging Software	7
Digital X-ray Sensors	8
Sensor Boot	8
Mouth Props	8
DC X-ray Generator	10
Handheld X-ray Generator	11
Complete Imaging System	11
CR Reader	12
DTX IO Camera	13

DENTAL STATIONS

Mobile	16 - 19
Wall Mount	20 - 22
Table Mount	22 - 23
Countertop	24 - 25
Compressors	28 - 29
Water Systems	30 - 31
Vacuum System	31

HANDPIECES

High-speed	34 - 39
Low-speed	40 - 42
Parts and Accessories	42 - 43

SCALERS AND MICROMOTORS

Scalers	46 - 49
Cavitron Inserts	49
Scaler Tips	50
Micromotors	51

CONSUMABLES

Prophy Paste	54
Prophy Cups	54
Fluoride	54
Prophy Angles	55 - 56
Contra Heads	57

HAND INSTRUMENTS

Scalers	70
Curettes	71
Excavators	72
Elevators	72 - 73
Winged Elevators	74 - 76
Periosteal Elevators	77
Luxators	78
Levitators	78
Scissors	79
Mouth Mirrors	79
Probes and Explorers	79 - 80
Forceps	80 - 81
Needle Holders	82
Root Tip Picks	82
Instrument Trays	82
Mouth Gags	83
Retractors	83 - 85
Rodent Instruments	86 - 87
Instrument Kits	88 - 93

ADVANCED DENTAL CARE

Bonded Sealant Kit	96
Curing Light	97
Restorative	98 - 104
Pluggers	105
K-reamers	106
H-files	106
Gutta Percha Points	107
Absorbant Paper Points	108
Sharpening Instruments	108
Impression Products	109

EDUCATIONAL RESOURCES

Mouth Models	112
Dental Charts	113

Click on item or category
to quickly view section.

DENTAL IMAGING

DTX Imaging System
DTX Imaging Software
Digital X-ray Sensors
Sensor Boot
Mouth Props
DC X-ray Generator
Handheld X-ray Generator
Complete Imaging System
CR Reader
DTX IO Camera



DC GENERATOR

DIGITAL SENSORS

HANDHELD GENERATOR

DTX IMAGING SOFTWARE

DTX Imaging System

Designed specifically for veterinary dentistry, the [DentalAire DTX Imaging System](#) includes essential products needed to capture high-quality dental images of patients.

With more than half of a patient's tooth below the gumline, taking dental x-rays with digital sensors enhances patient care and helps clinics properly diagnose and treat each patient.

DentalAire's DTX Imaging System is a bundled imaging package that was created to enhance your imaging needs, with lifetime phone and web-based support.

- Different size options for digital sensors
- Mouth prop kit to prevent accidental bites to the sensor
- Petsoft™ Pro for easy retraction and sensor placement
- Sensor boots for size 2 sensor to use as a positioning aid (only included with size 2 sensor)
- Sensor sleeves to protect sensors
- Innovative and customizable DTX Dental Software



Item Number	Description
DTP7000	DTX Digital Imaging System with Size 1 and 2 sensor
DTP70001	DTX Digital Imaging System with Size 1 sensor
DTP70002	DTX Digital Imaging System with Size 2 sensor
DTP700022	DTX Digital Imaging System with Size 4c sensor

DTX Imaging Software

- Direct e-mailer, image burner, and auto-route to multiple PACS
- Go from examination image to extraction procedure
- Auto full screen after capture
- Anatomically correct display on capture and view screen
- Advanced image comparison and image tagging functionality
- Efficient and improved installation process
- Real-time filters, advanced image labeling and enhancements
- Unique customization features
- Unlimited phone and web-based support
- Unmatched DICOM capabilities, including MWL query and auto-routing to PACS to quickly access and customize patient reports

DTX Software




Image Comparison



Item Number	Description
DTP700012	DTX Software License
DTP70004	DCV Software for DTX (software only)

Contact us
to schedule a live demo.

Patient Report




DentalAire
17742 Mitchell North, Suite B
Irvine, CA 92614
(949) 866-0881


Patient: JOHNSON, FLUFFY
Date: June 15, 2015



Right upper side showing severe rotation of the left incisors around the teeth.



The back of the mouth, showing severely rotated teeth in the back of the mouth.



View of the two lower canine teeth, showing severe loss of bone.



Extrusion of the left lower side, showing destruction of the teeth and bone.

Fluffy presented for a dental cleaning and radiographs of all her teeth. She had been scheduled to undergo an extraction for a loose tooth, which showed no signs of bone. The tooth was removed, and the extraction site was sutured. She was given a prescription for pain medication and instructed to return for a follow-up visit. The tooth was removed, and the extraction site was sutured. She was given a prescription for pain medication and instructed to return for a follow-up visit. The tooth was removed, and the extraction site was sutured. She was given a prescription for pain medication and instructed to return for a follow-up visit.

Please call if you have any questions or concerns.

DTX Digital Sensors

Available in 3 different sizes, the **DTX Digital Sensors** help produce high-quality images while also improving the overall quality of care for each patient.

- Direct USB eliminates connections and improves durability
- Waterproof and hermetically sealed housing is cold sterilizable and eliminates damage from bodily fluids
- High-quality CMOS chip, fiber optic plate, and CSI scintillator for superior image quality
- Patented top corner shape for easy positioning on molars
- Patented photodiode to automatically control over exposure
- DentalAire Sensor Replacement Program
- 3 year warranty*



DTP70008



DTP70009



DTP700020

Item Number	Description
DTP70008	DTX Size 1 Sensor
DTP70009	DTX Size 2 Sensor
DTP700020	DTX Size 4c Sensor

DentalAire DTX Sensor Replacement Program

- Universal eligibility, open to all clinics, regardless of sensor brand
- If the existing sensor is not covered by warranty, we offer a **NEW** replacement sensor at a discounted price, with no expiration date.
- If your sensor malfunctions due to manufacturers' workmanship during the 3-year warranty period, we will replace the sensor under our standard warranty replacement process.
- Sensor is latest model (not refurbished) with new warranty.
- No expiration date
- Discounted price from suggested retail with trade-in
- Sensor Replacement Program covers: accidental damage to sensor, sensor casing or cable, unusual wear and tear, improper cleaning or disinfecting of the sensor, cable, or cable connector


Item Number	Description
DTP700010	Size 1 Replacement Sensor
DTP700011	Size 2 Replacement Sensor
DTP700021	Size 4c Replacement Sensor

Sensor Boot

Designed to help protect your size 2 sensor and aid in your positioning. Autoclavable and sterilizable.

- Bright blue color for placement detection
- Grip for traction between teeth

Item Number	Description
DTP61613	DTX size 2 sensor boot



**digital sensor not included*

Sensor Sleeves

- Sleeve acts as a barrier between the sensor and patient's mouth
- Helps prevent scratches to the sensor
- Disposable and easy application



DTP61614

Item Number	Description
DTP61610	Size 2 Sensor Sleeves, 500 per box
DTP61614	Size 4c Sensor Sleeves, 500 per box

Mouth Props

- Used during a digital dental imaging procedure
- Provides bite prevention
- Autoclavable and sterilizable
- Available in sizes small, medium, and large

Item Number	Description
DTP61615	Small Mouth Prop
DTP61616	Medium Mouth Prop
DTP61617	Large Mouth Prop
DTP61618	Mouth Prop Kit, Set of Three (small, medium, and large)



DC X-ray Generator

The **DentalAire DC X-ray Generator's** ergonomic and modern design was built for veterinary demand. With high-resolution imaging, the DC X-ray Generator enhances patient care and provide clear images.

As clinic space and design can vary, the DC Generator is available in both a wall-mount and mobile version to meet your exact needs.

- Shorter exposure times and reduced dose
- Veterinary-specific selections including exotic and extremity
- 0.4 mm focal spot
- Full-size LCD display with adjustable brightness
- Adjustable kVp for maximum contrast adjustment
- High-frequency DC technology
- Clear tube cone allows for grasp and detection of sensor
- Long reach arm with frictionless joints
- Works seamlessly with DentalAire sensors
- Exposure can be made from the touch screen or wired exposure button
- Optional laptop tray (sold separately)
- 2 year warranty

Mobile

Item Number	Description
DTP50001	Mobile DC Generator
DTP61305	Mobile x-ray on-site installation
DTP50006	Laptop Tray for Mobile DC Generator

Wall Mount

Item Number	Description
DTP50002	DC Generator with 30 cm extension arm
DTP50003	DC Generator with 60 cm extension arm
DTP50004	DC Generator with 80 cm extension arm
DTP50005	DC Generator with 100 cm extension arm
DTP61304	Wall mount x-ray on-site installation
DTP60808	Laptop Tray for Wall Mount DC Generator

Dental Positioning Training

Item Number	Description
DTP700018	On-site Dental Positioning Training (4 hours)



Handheld Portable X-ray Generator

The **DentalAire Handheld Portable X-ray Generator's** lightweight design increases mobility and provides high-quality images. The included neck strap helps the operator when using the generator to prevent falls and accidental damage, while also reducing arm fatigue.

- 0.4 mm focal spot
- Ergonomic camera-type design
- 4.4 pounds
- Permanently attached backscatter shield
- Included storage case
- 1 year warranty*



Item Number	Description
DTP61307	Handheld X-ray Generator

Complete Imaging Systems

Complete your imaging system with the **DentalAire Complete Imaging System**. It is designed to seamlessly work together in your clinic to improve your workflow and optimize performance. DentalAire also offers technical support for all components of the system.

Included in the Complete Imaging System is:

- X-ray generator of your choice
- DTX software
- Size 2 sensor
- Laptop
- Petsoft™ Pro, Mouth Props , Sensor Boot
- Unlimited phone and web-based support



Wall Mount

Item Number	Description
DTP60811	30 cm arm wall mount generator, laptop, size 2 sensor, DTX software
DTP60812	60 cm arm wall mount generator, laptop, size 2 sensor, DTX software
DTP60813	80 cm arm wall mount generator, laptop, size 2 sensor, DTX software
DTP60816	100 cm arm wall mount generator, laptop, size 2 sensor, DTX software

Mobile

Item Number	Description
DTP61310	Mobile DC Generator, laptop, size 2 sensor, DTX software

Handheld

Item Number	Description
DTP60815	Handheld Generator, laptop, size 2 sensor, DTX software

DTX CR Reader

The **DentalAire DTX CR Reader** helps optimize the workflow of your practice with rapid, high-quality readings of reusable dental imaging plates.

- Easy-to-use LCD touch screen interface
- Imaging plates included: (1) - Size 4c and (4) - Size 2
- Included DTX software
- Easy-to-position bite-wing & intraoral imaging plates
- Built-in erase function
- Software optimizes acquisition, processing & management
- Supports both Ethernet & USB connections
- Auto-calibration ensures stability & consistency



Item Number	Description
DTP710026	DentalAire CR Reader with DTX software

Imaging Plates

Item Number	Description
DTP710027	Size 0 Imaging Plate Kit Contains: 4 IP
DTP710029	Size 1 Imaging Plate Kit Contains: 4 IP
DTP710022	Size 2 Imaging Plate Kit Contains: 4 IP
DTP710023	Size 4c Imaging Plate Kit Contains: 4 IP



Hygienic Bags

- Help protect the imaging plates
- Easy to use - loading, sealing, and opening

Item Number	Description
DTP710028	Hygienic Bags Size 0 box of 100 pcs
DTP710030	Hygienic Bags Size 1 box of 100 pcs
DTP710024	Hygienic Bags Size 2 box of 100 pcs
DTP710025	Hygienic Bags Size 4c box of 100 pcs



DTX IO Camera

The **DTX IO Camera** has a professional curved optical lens and a curved design for easier photo shots. It captures crisp, clear images while eliminating distortion.

- Interfaces with DTX Dental Imaging System
- 6 LED lights for high-illumination
- USB 2.0 with 2.5m cable length
- 1 year warranty

Item Number	Description
DTP00191	DTX IO Camera



DTX Laptop

Item Number	Description
LAPTOP	DentalAire DTX Laptop



Sani-Cloth® Wipes

- Free from alcohol, fragrance, bleach, acid, phenol, acetone, and ammonia
- Effective against 45 clinically relevant microorganisms
- Compatible with a broad range of hard non-porous surfaces and equipment
- Wipe size: 6" x 6.75"
- 160 wipes per can



Item Number	Description
DTP10487	Sani-cloth ® AF3 Germicide Wipes

DENTAL STATIONS

Mobile
Wall Mount
Table Mount
Countertop
Compressors
Water Systems
Vaccuum System



PRESTIGE

ULTIMATE

COUNTERTOP III

DIPLOMATE STANDARD

Prestige Dental Station

The **Prestige Dental Station** is a popular mobile dental station that is known for its reliability, internal storage, and pneumatic height adjustment.

- Equipped with high-quality handpieces to enhance patient care:
New LED high-speed swivel handpiece, low-speed handpiece, LED Piezo Scaler, and air/water syringe
- All dental stations come with the **Dental Vet Package** that includes mint prophylaxis paste, PetSoft Pro™ snap-type prophylaxis cups, snap-on prophylaxis angles, bur starter pack, o-ring for tank drain, and handpiece cleaner and lubricant.
- Offers telescopic height adjustment from an operating height of 38 1/4" down to 30" for storage purposes when not in use
- Hassle-free quick disconnect pressurized water connection
- Standard Midwest 4-hole connection for optimal versatility, performance, and convenience
- Oil, oil-less, or no compressor option designed to meet a clinic's specific needs and requirements
- Asepsis design with controls located on the side panel
- Dimensions: 38 1/4" - 30" h x 22 1/2" w x 24 1/4" - 22" d
- 2 year warranty*



DTP00511

Item Number	Description
DTP00511	Prestige Dental Station with Oil Compressor
DTP00511-OL	Prestige Dental Station with Oil-less Compressor
DTP00511-2HP	Prestige Dental Station with Two LED High-speed Handpieces
DTP00512	Prestige Dental Station with No Compressor



Ultimate Dental Station

The **Ultimate Dental Station** has a large workspace area that includes a removable stainless-steel tray which can be used to place dental instruments and supplies.

- Equipped with high-quality handpieces to enhance patient care: **New** LED high-speed swivel handpiece, low-speed handpiece, Piezo Scaler, and air/water syringe
- All dental stations come with the **Dental Vet Package** that includes mint prophy paste, PetSoft Pro™, snap-type prophy cups, snap-on prophy angles, bur starter pack, o-ring for tank drain, and handpiece cleaner and lubricant.
- Offers manual height adjustment from an operating height of 40" down to 26.5" for storage purposes when not in use
- Hassle-free quick disconnect pressurized water connection
- Standard Midwest 4-hole connection for optimal versatility, performance, and convenience
- Oil, oil-less, or no compressor option designed to meet a clinic's specific needs and requirements
- Dimensions: 40" - 26.5" h x 22" w x 22" d
- 1 year warranty



Item Number	Description
DTP00509	Ultimate Dental Station with Oil Compressor
DTP00509-OL	Ultimate Dental Station with Oil-less Compressor
DTP00510	Ultimate Dental Station with No Compressor



Basic Dental Station

The **Basic Station** can easily move around your veterinary practice when space is limited. With this station, a clinic can use their air source or it can be paired with the DentalAire compressor.

Its portability and simplicity make it ideal for clinics looking to maximize small spaces.

- Equipped with high-quality handpieces to enhance patient care: **New** LED high-speed swivel handpiece, low-speed handpiece, built-in LED ultrasonic scaler, and air/water syringe
- All dental stations come with the **Dental Vet Package** that includes mint prophylaxis paste, PetSoft Pro™, snap-type prophylaxis cups, snap-on prophylaxis angles, bur starter pack, o-ring for tank drain, and handpiece cleaner and lubricant.
- One compact 2-liter water bottle system to maximize convenience and portability
- Hassle-free quick disconnect pressurized water connection
- Standard Midwest 4-hole connection for optimal versatility, performance, and convenience
- Variable speed foot control to improve ergonomics and help enhance performance with procedures and dental cleanings
- 1 year warranty




Item Number	Description
DTP00497	Basic Dental Station

Dentaire II Station

The **Dentaire II Station** includes a removable stainless-steel tray that can be used to place dental instruments and supplies.

- Equipped with high-quality handpieces to enhance patient care: high-speed handpiece, low-speed handpiece, and air/water syringe
- All dental stations come with the **Dental Vet Package** that includes mint prophy paste, PetSoft Pro™, snap-type prophy cups, snap-on prophy angles, bur starter pack, o-ring for tank drain, and handpiece cleaner and lubricant.
- Variable speed foot control to improve ergonomics and help enhance performance with procedures and dental cleanings
- Standard Midwest 4-hole connection for optimal versatility, performance, and convenience
- Oil or no compressor options designed to meet a clinic's specific needs and requirements
- 1 year warranty



DTP00495

Item Number	Description
DTP00495	Dentaire II Station with Oil Compressor
DTP00494	Dentaire II Station with No Compressor

Diplomate Standard Wall Mount

The **Diplomate Wall Mount Dental Station** is our most advanced wall-mount dental system. With a flip-top delivery system, it allows for the internal components to be easily accessible. The flex arm system provides more flexibility and efficiency when in need of a space-saving solution.

- Equipped with high-quality handpieces to enhance patient care:
New LED high-speed swivel handpiece, low-speed handpiece, LED Piezo Scaler, and air/water syringe
- All dental stations come with the **Dental Vet Package** that includes mint prophylaxis paste, PetSoft Pro™, snap-type prophylaxis cups, snap-on prophylaxis angles, bur starter pack, o-ring for tank drain, and handpiece cleaner and lubricant.
- Side-mounted control adjustments
- 3-way syringe with quick change autoclavable tip
- Asepsis flex arm system with pneumatic brake, 52" reach, 27" vertical travel
- Soft touch variable speed foot control
- Dual 2-liter water bottles with wall mount or post mount option
- Nitrogen or compressor air-driven
- 2 year warranty*



DTPJ-3300

Item Number	Description
DTPJ-3300	Diplomate Wall Mount Standard Style with Oil Compressor
DTPJ-3300-OL	Diplomate Wall Mount Standard Style with Oil-less Compressor
DTPJ-3300-NC	Diplomate Wall Mount Standard Style with No Compressor
DTPJ-3300-P	Diplomate Wall Mount Standard Style with Plumbed Air/Water Regulators
DTPJ-3300-PC	Diplomate Wall Mount Standard Style with Plumbed Air/Water Regulators Plus Compressors



Diplomate Basic

The **Diplomate Basic** is a wall-mount dental station that provides flexibility and efficiency for veterinary practices in need of space-saving solutions. With customizable scissor arms, it can efficiently fit into your clinic.

- Equipped with high-quality handpieces to enhance patient care: **New** LED high-speed swivel handpiece, low-speed handpiece, LED Piezo Scaler, and air/water syringe
- All dental stations come with the **Dental Vet Package** that includes mint prophylaxis paste, PetSoft Pro™, snap-type prophylaxis cups, snap-on prophylaxis angles, bur starter pack, o-ring for tank drain, and handpiece cleaner and lubricant.
- Variable speed foot control to improve ergonomics and help enhance performance with procedures and dental cleanings
- Standard Midwest 4-hole connection
- Included Vet Package
- One (1) - 2-liter water bottle
- 1 year warranty



Item Number	Description
DTP00705	Diplomate Wall Mount Basic with Oil Compressor

Three Handpiece Wall Mount Unit

- **New** LED swivel high-speed handpiece, low-speed handpiece, Piezo Scaler, and air/water syringe
- Variable speed foot control
- Mix and match any two arm lengths (6", 12", 16") or the standard option that includes 2 - 12" arms
- Available in scissor configuration or single-arm
- Included Vet Package
- One (1) - 2-liter water bottle
- 1 year warranty



Item Number	Description
DTP00603	Three Handpiece Wall Mount Unit with Oil Compressor
DTP00603-OL	Three Handpiece Wall Mount Unit with Oil-less Compressor
DTP00603-LED	Three Handpiece Wall Mount Unit with LED Scaler
DTP00603-LED-NC	Three Handpiece Wall Mount Unit with LED Scaler and No Compressor
DTP00603-NHP-NC	Three Handpiece Wall Mount Unit with No Handpieces and No Compressor

Two Handpiece Wall Mount System

- Equipped with high-quality handpieces to enhance patient care: high-speed push button handpiece, low-speed handpiece, and air/water syringe
- One compact 2-liter water bottle system to maximize convenience and portability
- Standard Midwest 4-hole connection for optimal versatility and performance
- Variable speed foot control to improve ergonomics and help enhance performance with procedures and dental cleanings
- Included Vet Package
- 1 year warranty



Item Number	Description
DTP00600	Wall-Mounted Two Handpiece System without Compressor
DTP00601	Wall-Mounted Two Handpiece System with Compressor

Two Handpiece Table Mount System

The **DentalAire Table Mount System** can be clamped onto a table or other similar surfaces. Because this unit can clamp onto things such as tables, it is a popular style for labs.

- Stainless steel tray
- High-speed push button handpiece, low-speed handpiece, and water/air syringe
- One (1) - 750ml water bottle
- Variable speed foot control
- Included Vet Package
- 1 year warranty



Item Number	Description
DTP00607	Table Mounted Two Handpiece System without Compressor
DTP00608	Table Mounted Two Handpiece System with Compressor

DentalAire Olympic Wet Vet Table Dental Station

The **DentalAire Olympic Wet Table Dental Machine** is specifically made for the Olympic Wet Table. It easily attaches to the wet table to aid in dental cleanings.

- Equipped with high-quality handpieces to enhance patient care: **New** LED swivel high-speed handpiece, low-speed handpiece, LED Piezo scaler and air/water syringe
- Variable speed foot control to improve ergonomics and help enhance performance with procedures and dental cleanings
- Two (2) - 2 liter water bottles
- Standard Midwest 4-hole connection for optimal versatility, performance, and convenience
- Included Vet Package
- Oil compressor



DTP00700



*Olympic wet table not included

Item Number	Description
DTP00700	Olympic Wet Table Dental Machine
DTP00700-NC	Olympic Wet Table Dental Machine with No Compressor



Countertop II Station

The **DentalAire II Countertop Unit** is a two-handpiece system.

- Equipped with high-quality handpieces to enhance patient care: high-speed push button handpiece, low-speed handpiece, and air/water syringe
- All dental stations come with the **Dental Vet Package** that includes mint prophylaxis paste, PetSoft Pro™, snap-type prophylaxis cups, snap-on prophylaxis angles, bur starter pack, o-ring for tank drain, and handpiece cleaner and lubricant.
- One compact 750ml water bottle system to maximize convenience and portability
- Durable handle to improve portability and make it easier to carry dental station
- Standard Midwest 4-hole connection for optimal versatility and performance
- Variable speed foot control to improve ergonomics and help enhance performance with procedures and dental cleanings
- 1 year warranty



Item Number	Description
DTP00500	DentalAire II Countertop Station with Oil Compressor
DTP00500-OL	DentalAire II Station with Oil-less Compressor

Countertop III Station

The **DentalAire III Station** is a three-handpiece system.

- Equipped with high-quality handpieces to enhance patient care: high-speed handpiece, low-speed handpiece, Piezo Scaler, and air/water syringe
- Variable speed foot control to improve ergonomics and help enhance performance with procedures and dental cleanings
- Standard Midwest 4-hole connection for optimal versatility and performance
- One compact 750ml water bottle system to maximize convenience and portability
- Durable handle to improve portability and make it easier to carry dental station with oil compressor
- Dimensions: 19" h x 11" w x 18" l
- Included Vet Package
- 1 year warranty



Item Number	Description
DTP00504	DentalAire III Countertop Station with Push Button High-speed Handpiece
DTP00504-LED	DentalAire III Station with LED Swivel High-speed Handpiece and LED Scaler





DTPJ-3300



DTP00511



DTP00705



DTP00504

DENTAL STATION COMPARISON CHART

Mobile Dental Stations

Item Number	Name	Scaler	High-speed	Low-speed	Water Bottle	Compressor
DTP00494	Mobile Cart II	No	Push Button	Yes	2 Liters	No
DTP00495	Mobile Cart II	No	Push Button	Yes	2 Liters	Oil
DTP00497	Basic	Built in LED Ultrasonic	LED Swivel PB	Yes	2 Liters	No
DTP00509	Ultimate	LM Piezo	LED Swivel PB	Yes	2 - 2 Liters	Oil
DTP00509-OL	Ultimate	LM Piezo	LED Swivel PB	Yes	2 - 2 Liters	Oil-less
DTP00511	Prestige	LED Piezo	LED Swivel PB	Yes	2 - 2 Liters	Oil
DTP00511-OL	Prestige	LED Piezo	LED Swivel PB	Yes	2 - 2 Liters	Oil-less
DTP00511-2HP	Prestige	LED Piezo	LED Swivel PB	Yes	2 - 2 Liters	Oil
DTP00512	Prestige	LED Piezo	LED Swivel PB	Yes	2 - 2 Liters	No

Countertop Dental Stations

Item Number	Name	Scaler	High-speed	Low-speed	Water Bottle	Compressor
DTP00500	Countertop II	No	Push Button	Yes	750 ml	Oil
DTP00500-OL	Countertop II	No	Push Button	Yes	750 ml	Oil-less
DTP00504	Countertop III	LM Piezo	Push Button	Yes	750 ml	Oil
DTP00504-LED	Countertop III	LED Piezo	LED Swivel PB	Yes	750 ml	Oil

Wall Mounted Dental Stations

Item Number	Name	Scaler	High-speed	Low-speed	Water Bottle	Compressor
DTPJ-3300	Diplomate Standard	LED Piezo Ultrasonic	LED Swivel PB	Yes	2 - 2 Liters	Oil
DTPJ-3300-OL	Diplomate Standard	LED Piezo Ultrasonic	LED Swivel PB	Yes	2 - 2 Liters	Oil-less
DTPJ-3300-NC	Diplomate Standard	LED Piezo Ultrasonic	LED Swivel PB	Yes	2 - 2 Liters	No
DTP00603	3 Handpiece	LM Piezo	LED Swivel PB	Yes	2 Liters	Oil
DTP00603-OL	3 Handpiece	LM Piezo	LED Swivel PB	Yes	2 Liters	Oil-less
DTP00603-LED	3 Handpiece	Built in LED	LED Swivel PB	Yes	2 Liters	Oil
DTP00603-LED-NC	3 Handpiece	Built in LED	LED Swivel PB	Yes	2 Liters	No
DTP00603-NHP-NC	3 Handpiece	No	No	No	2 Liters	No
DTP00600	Wall Mount II	No	Push Button	Yes	2 Liters	No
DTP00601	Wall Mount II	No	Push Button	Yes	2 Liters	Oil
DTP00705	Diplomate Basic	LED Piezo Ultrasonic	LED Swivel PB	Yes	2 Liters	Oil

Table Mounted Dental Stations

Item Number	Name	Scaler	High-speed	Low-speed	Water Bottle	Compressor
DTP00607	Table Mount II	No	Push Button	Yes	750 ml	No
DTP00608	Table Mount II	No	Push Button	Yes	750 ml	Oil
DTP00700	Olympic Wet Vet...	Built in LED	LED Swivel PB	Yes	2 - 2 Liters	Oil
DTP00700-NC	Olympic Wet Vet...	Built in LED	LED Swivel PB	Yes	2 - 2 Liters	No

*All equipment above apply the following features: High-Speed, Low-Speed, Water & Air, unless stated otherwise.

Oil Air Compressor 110V

- H.P. 1/2
- Voltage: 110/60, 230/50-60
- Output: 2 CFM
- Max Pressure: 114 PSI
- Operating Pressure: 84-114 PSI
- Noise Level: 40 Decibels
- Tank Size: 2.25 Gallons
- Dimensions: 17x7.9x17.7
- Weight: 58 lbs



Item Number	Description
DA50-9D	DentalAire Oil Compressor

Oil Compressor with Fans

- H.P. 2 X 1/2
- Voltage: 110/60, 230/50-60
- Output: 4.15 CFM
- Max Pressure: 114 PSI
- Operating Pressure: 84-114 PSI
- Noise Level: 42 Decibels
- Tank Size: 6 Gallons
- Dimensions: 25.6x14.2x21.6
- Weight: 100 lbs



Item Number	Description
DA100-24	1 HP Sil-Air Compressor with Fans

Oil-less Compressor

- Noise level: 56 Decibels
- Overall dimensions 17.5 x 9 x 15.5
- Tank .93 Gal
- CFM: 3.25
- 115 V/50 - 60 Hz or 220 V/50 - 60 Hz



Item Number	Description
DA50-9D3	DentalAire Oil-less Compressor

Oil-less Compressor with Dryer and Silencing Cabinet

- Noise Level: 52 Decibels
- Tank Size: 13 Gallons
- Dimensions: 29x23x32
- Weight: 207 lbs

Item Number	Description
DA2/50/379	1.5 HP Oil-less Compressor with Dryer and Silencing Cabinet 110V



Oil-less Compressor with Dryer

- H.P: 1.5
- Voltage: 110/60, 230/50, 230/60
- Output: 6.7 CFM / 190 LPM
- Max Pressure: 120 PSI
- Operating Pressure: 90-120 PSI
- Noise Level: 72 decibels
- Tank Size: 13 gallons
- Dimensions: 18 x 26 x 35
- Weight: 119 lbs

Item Number	Description
DA2/50/37	1.5 HP Oil-less Compressor with Dryer 110V



Oil Compressor Maintenance Kits

The **Compressor Maintenance Kits** are ideal for all DentalAire oil-compressors to help maintain and ensure the compressor is running efficiently.

- C00PM and C00PM-2: Contain 22 oz oil bottle, drain tank o-ring, metal glass o-ring, oil fill plug o-ring, air intake filter element
- C0164-SMCKIT: Filter Element/Filter Regulator, Coalescent Replacement Filter, Pressure Bowl O-ring SMC

Item Number	Description
C00PM	Compressor Maintenance Kit
C00PM-2	Compressor Maintenance Kit with SMC Filter Kit
DTP11310	Compressor Oil 22 oz Only
C0164-SMCKIT	SMC Filter Kit



C00PM

Water Distiller

A reliable source for high-quality distilled water.

- Voltage: 100 - 240V
- Distiller water: 1.5 L/H
- Frequency: 50 - 60 HZ
- Weight: 3.5kgs (7.72 lbs)
- 1 year warranty



Item Number	Description
DTP00650	Water Distiller 4L

Pump Up Water Bottle System

Portable water tank with built-in hand pump for distilled water or medicated solution.

- 1-gallon capacity

Item Number	Description
DTP00773	Pump Up Water Bottle System



Water Bottle Cartridge

Dental unit waterline treatment cartridges that deliver clean water to staff for 1 year or 240L of water.

- Elemental iodine in DentaPure contains no allergenic proteins
- Safe and effective uninterrupted waterline treatment



DTP365M

Item Number	Description
DTP365B	365 Independent Water Bottle Cartridge
DTP365M	365 Municipal (City) Water Bottle Cartridge



DTP365B

DENTAL STATIONS | WATER AND VACUUM

Water Bottles

Item Number	Description
A10-002	Water Bottle 750ml
A110-003	Water Bottle 2 Liters



Sterident Vacuum System

Efficient, adjustable motor speed can evacuate indefinitely at the rate of 75 cfm or can be adjusted.

- Compact (13" x 16" x 26")
- Easy maintenance
- Non-corrosive canister
- Cleans efficiently
- Equipped with automatic overflow shut-off
- Motor: 1 year warranty

Item Number	Description
DA8009PD	Sterident Port-Dent II Mobile Evacuation Unit



AVS Vacuum System

- Self-contained air-powered vacuum system
- For optimal use, dual user compressor should be ordered

Item Number	Description
AVS	AVS Vacuum System



HANDPIECES



LED HIGH-SPEED

FIBER-OPTIC HIGH-SPEED

LOW-SPEED

HANDPIECE LUBRICANT

LED Swivel High-speed Handpiece

- Triple water spray cools the entire surface of the bur to avoid overheating
- Built-in generator
- Push button
- Innovative micro-generator
- 410,000 RPM

Item Number	Description
DA400LEDMH	LED Swivel High-speed Handpiece, Mini Head

Parts and Accessories

Item Number	Description
DA400MHTURB	Turbine for Swivel LED Highspeed Mini Head
DA400LEDMH-BC	Mini Head Back Cap
DA400LEDMH-SC	Swivel Coupler for MH LED
DA400LEDMH-ORS	Swivel Coupler O-ring Set
DA400LEDMH-HDCT	Heavy Duty Cap Tool for MH Swivel LED Drill
DTPLED-CPT	Standard Cap Tool
DA400SP	Shipping Post



LED High-speed Handpiece

- Triple water spray cools the entire surface of the bur to avoid overheating
- Built-in generator
- Push button
- Innovative micro-generator
- 410,000 RPM

Item Number	Description
DA400LED	LED High-speed Handpiece, Standard Head

Parts and Accessories

Item Number	Description
DA400LED LG TURB	Turbine for LED Highspeed
DA400LED CP	Push Button Cap
DA400LED-HDCT	Heavy Duty Cap Tool for LED Drill
DA400LED-CPT	Standard Cap Tool
DA400SP	Shipping Post



Fiber-optic Swivel High-speed Handpiece

- Fiber-optic light
- Push button
- 400,000 RPM

Item Number	Description
-------------	-------------

DA400AFOPB-SW	Fiber-optic Swivel High-speed Handpiece
---------------	---

Parts and Accessories

Item Number	Description
-------------	-------------

DANSK-A485	Turbine
DA400RCA	Turbine and Push Button Cap
VKS-5	Kavo Coupler for Glass Rod
DA400SORK	Swivel HP Kavo Coupler O-ring Kit
DA100-CAPTOOL	Cap Tool
DA400SP	Shipping Post
DANSK-A100	Push Button Cap



Fiber-optic High-speed Handpiece

- Fiber-optic light
- Push button
- 400,000 RPM

Item Number	Description
DA400AFOPB	Fiber-optic High-speed Handpiece

Parts and Accessories

Item Number	Description
DANSK-A485	Turbine
DA400RCA	Turbine and Push Button Cap
DANSK-A100	Push Button Cap
DA100-CAPTOOL	Cap Tool
DA400SP	Shipping Post



Push Button High-speed Handpiece

- 400,000 RPM
- Push button
- Single water spray

Item Number	Description
DA400APB	High-speed Push Button

Parts and Accessories

Item Number	Description
DANSK-A485	Turbine
DA400RCA	Turbine and Push Button Cap
DANSK-A100	Push Button Cap
DA100-CAPTOOL	Cap Tool
DA400SP	Shipping Post



Swivel Push Button High-speed Handpiece

- 400,000 RPM
- Push button
- Single water spray

Item Number	Description
DA400SWPB	High-speed Swivel Push Button

Parts and Accessories

Item Number	Description
DAMCT-A445	Turbine
MC	Dentex Non-Fiber Optic Coupler
DA400SWPBCAP	Dentex M3MB-SLS P.B. Cap
DA100-CAPTOOL	Cap Tool
DA400SP	Shipping Post

Standard High-speed Handpiece

- 400,000 RPM
- Single water spray
- Manual chuck

Item Number	Description
DA400A	Standard High-speed Handpiece

Parts and Accessories

Item Number	Description
DA400RC	Turbine
DA400CP	Standard Replacement Cap
DA100-CAPTOOL	Cap Tool
DA400SP	Shipping Post
DA400CH	Replacement Chuck for Standard Turbine



Standard Fiber-Optic High-speed Handpiece

- 400,000 RPM
- Single water spray
- Manual chuck

Item Number	Description
DA400AFO	Standard Fiber-Optic High-speed Handpiece

Parts and Accessories

Item Number	Description
DA400RC	Standard Turbine
DA400CP	Standard Replacement Cap
DA100-CAPTOOL	Cap Tool
DA400SP	Shipping Post

Hygienist Low-speed

- 5,000 RPM
- U-type nose
- ISO 4-holes
- Push/Pull chuck system



Item Number	Description
DA600HA	Hygienist Low-speed Handpiece

Lube-Free Low-speed

- Speed-adjusting ring controls the speed and direction of the motor
- Swivels 360 degrees
- Straight nose cone
- 5,000 RPM



Item Number	Description
DA700A	Lube-Free Low-speed Handpiece

Parts and Accessories

Item Number	Description
DAMLS137A	Slow-speed Swivel Tip for DA700A

Standard Low-speed

- Opening and closing chuck system and reverse
- ISO 4-holes
- U-type nose
- 5,000 RPM



Item Number	Description
DA600A	Standard Low-speed Handpiece

Low-speed

- Accepts standard straight handpiece burs, contra-angles, and disposable prophyl angles
- 360° swivel feature prevents tube twisting and back-end drag
- Forward and reverse rotational control
- ISO 4-holes
- 5,500 RPM



Item Number	Description
DA500A	Low-speed Handpiece

E-Type Straight Nose Cone

- 40,000 RPM
- Applicable bur: HP type ISO standard

Item Number	Description
DA2500SNC	E-type Straight Nose Cone



E-type Air Motor

- 20,000 RPM
- Applicable bur: HP type ISO standard

Item Number	Description
DA2500EM	E-type Air Motor



Contra Angle - Latch Type

- 15,000 RPM
- U-type contra angle
- Applicable bur: CA type ISO standard

Item Number	Description
DA553UL	Contra Angle - Latch Type



Contra Angle - Friction Grip

- 15,000 RPM
- U-type coupling
- Applicable bur: FG type ISO standard

Item Number	Description
DA557FG	Contra Angle - Friction Grip



E-type Contra Angle - Latch Type

- 20,000 RPM
- ISO standard (E-type)
- Applicable bur: CA type ISO standard



Item Number	Description
DA2500ECA	E-type Contra Angle - Latch Type

Sonic Air Scaler Handpiece

- Controlled water spray system
- Helps remove deposits, stains, and calculus
- Includes 3 scaler tips




Item Number	Description
DA7000AS	Sonic Air Scaler Handpiece

Handpiece Lubricant

- Promotes smoother operating sound
- Helps the handpiece perform efficiently
- Extends the longevity of handpiece

Item Number	Description
DA100PL	Pen Handpiece Lubricant 1/4 oz
DA100JP	All-in-one Handpiece Cleaner and Lubricant
DA100MC	Air Motor Cleaner



DA100JP



DA100PL



DA100MC

Gaskets for Handpieces

Item Number	Description
DA500GAS	Gasket for Fiber-optic Handpiece
DA400GAS	Gasket for 4-hole Handpieces



DA400GAS



DA500GAS

Smart Kit®

- Removes water and biofilm buildup
- Minimizes delays and maximizes clinic productivity
- Prevents initial residue/particle build-up
- Fits Midwest standard connection and e-type contra-angle handpieces



What's Included:

1



Smart Cleaner®

- Tool for clearing blockages in handpieces and quick coupling

2



Smart Connect®

- Adaptor for contra-angle handpiece (E-type) with internal water

3



Smart Cupper®

- Specialized tool for isolated blockages in multi-hole handpieces

Item Number	Description
DTP45344	Smart Kit®, Handpiece Water Line Cleaner

Fiber-optic Bulb

Item Number	Description
DA400FOB	Fiber-optic bulb



Bur Changer

Item Number	Description
DA310BT	Handpiece Bur Changer



SCALERS AND MICROMOTORS



NEWTRON P5 WITH LED

LM ERGO VET SCALER

PIEZOTOME® CUBE

COCOON SPRAY

LM Ergo Vet Scaler

The **ErgoVet Scaler** offers ultrasonic solutions not only for traditional scaling and periodontics, but also for endodontics. It is an effective piezoelectric scaler in one versatile and ergonomic appliance.

The device's LED lights, advanced electronics, quality, and high-durability LM-DuraGradeMAX tips enhance the execution of procedures. Ergonomically designed LM-ErgoGrip handpieces with soft silicone offer a comfortable, relaxed grip.

- Removes biofilm and calculus on teeth and other dental work where the ultrasonic vibration is beneficial
- Piezo-electric technology
- Dynamic foot pedal for control and power adjustment
- Autoclavable handpiece
- LM-DuraGradeMax super steel tips
- Includes: handpiece and tubing, 3 scaling tips, 3 LM-ErgoGrips, and 3 torque wrenches
- Optimized LED lights
- Automatic cleaning function
- 270mm L x 140mm W x 165mm H
- 1 year warranty

Item Number	Description
-------------	-------------

DTP00901	ErgoVet Scaler with Built in Water Bottle
DTP00900	ErgoVet Scaler with Water Line Connector



Recommended LM Replacement Scaler Tips

Item Number	Description	Item Number	Description
DTP00784	Universal Tip #39	DTP00791	Interproximal #37
DTP00785	Flat Tip #32	DTP00793	Catheter Tip
DTP00786	Special / Perio Tip #33	DTP00794	Scaling tip #31 (Removes Gross Calculus)
DTP00788	Endo File Holder Tip 90°	DTP00796	Scaling Tip #PE-41PS (Super Slim, Perio)
		DTPTIPPAK	Set of 5 LM tips, 31, 32, 33, 37, 39

P7 Scaler

- Footpedal and Handpiece
- Scaling Tips (5): T1, T2, T4, T5, T3 Perio
- E1 (and 4 files) Endodontic function
- Torque Wrench
- Endo Chuck & Wrench
- Power Adapter
- Waterline (water source not included)
- 1 year warranty



DTP00798

Item Number	Description
DTP00799	P7 LED Ultrasonic Scaler (Scaler/Perio Endo)
DTP00798	P7 Ultrasonic Scaler (Scaler/Perio Endo)

Newtron® P5 with LED

Newtron P5 preserves the teeth and the soft tissue with controlled ultrasonic vibrations, specific quality of the steel of the tips, and irrigation that reaches the end of the tip. Linear and regular vibrations make treatments more comfortable for the patient.

- Flat glass surface
- Watertight and hygienic
- Removable power knob
- Adjustable, autoclavable handpiece holder
- Connects to water supply via Quick-Connect
- Unit contains: Newtron slim handpiece, white LED ring, 3 tips (#1,10Z, #H3), blue dynamometric wrench (2), green dynamometric wrench (1), autoclavable plastic box
- 2 year manufacturer warranty
- 90 day scaler tip manufacturer warranty



Item Number	Description
DTP00797	Newtron® P5 with LED

Newtron® Booster

Newtron® Booster preserves the teeth and the soft tissue. The perfectly controlled ultrasonic vibrations allow precise and gentle treatments. With the specific quality of the steel of the tips, it preserves the anatomy and enamel.

- Automatic and continuous frequency adjustment
- Linear and regular vibrations
- Total irrigation control
- Color Coding System that associates each tip color ring with one of the 4 power ranges
- Includes: Newtron slim handpiece, 3 tips (#1, #2, #10P), Blue dynamometric wrench (3), Autoclavable plastic box
- 2 year manufacturer warranty
- 90 day scaler tip manufacturer warranty



Item Number	Description
DTP00795	Newtron® Booster

Cocoon Spray

The **Cocoon Spray** brings high-performance compact equipment to veterinary dentistry. This model is equipped with piezoelectric technology.

- Combination scaler and electric polisher
- Rapidity: mechanical elimination of deposits assisted by cavitation
- Instrument Control: generating perfectly-controlled linear vibrations
- Power: self-regulated frequency (27 to 32KHz) for greater impact than pneumatic scalers
- Scaler: Compatible with both diamond-tipped and tungsten carbide burs
- Polisher: Spray cooling to avoid tissue necrosis and to remove debris
- Includes: #1 Universal, #2 Flat, #10P Shallow Pockets, Latch Type Contra Angle, Straight Nose Cone, Prophy Angle, Prophy Cups, Tip Wrench
- Performance: 40,000 RPM
- 2 year manufacturer warranty
- 90 day scaler tip manufacturer warranty



Item Number	Description
DTP24484V-K	Cocoon Spray, Combination Scaler/Electric Polisher

Piezotome® Cube

Piloted by advanced Newtron® technology, the **Piezotome® Cube** drastically improves the experience of dental procedures, while also improving the quality of care to patients.

- Faster, atraumatic and more predictable bone treatments
- Alternation between high and low-intensity frequencies
- Exclusive Tip Designs: Irrigation flow reaches the end of the tip regardless of orientation
- Minimal bone loss for a less invasive procedure
- Microbubbles releasing oxygen limiting the bleeding
- Autoclavable handpiece holder with 2 positions protecting from accidental fall
- Includes Bone surgery™ Kit and LC Extraction II Tip Kit
- 2 year manufacturer warranty*



DTPF50108

Bone surgery™ Kit



DTPF87509

- Intended for bone block grafting, cutting, scraping, and remodeling bone structures
- Included tips: BS1S, BS2L, BS2R, BS4, BS5, BS6
- 90 day scaler tip manufacturer warranty

LC Extraction II Tip Kit



DTPF87546

- Intended for extraction of teeth, fractured roots, impacted teeth, and root section
- Included tips: LC1, LC2, LC2L, LC2R, NINJA
- 90 day scaler tip manufacturer warranty

Item Number	Description
DTPF50108	Piezotome® Cube
DTPF87509	Bone surgery™ Kit for Piezotome®
DTPF87546	LC Extraction II Tip Kit

Cavitron Inserts

Item Number	Description
DTP10879	Beavertail Cavitron Insert Fine 25K
DTP10880	Universal Insert Fine 25K
DTP10881	Cavitron Insert Fine 30K



DTP10879



DTP10880



DTP10881

Acteon Satalec

Item Number	Description
DTPF00065	Apical Surgery AS3D Tip
DTPF00245	Perio #1S Slim Tip
DTPF00246	#1 Universal Tip
DTPF00247	#2 Voluminous Calculus Tip
DTPF00248	#3 Tip
DTPF00253	10p Shallow Pockets Tip
DTPF00254	#10Z Subgingival Tip
DTPF00255	3 Probe Urolithiasis Kit
DTPF00256	Urolithiasis Probe For Dogs
DTPF00257	Urolithiasis Probe For Cats
DTPF00366	H1 Curette Tip
DTPF00369	H3 Universal Currette Tip

Black Air Scaler

Item Number	Description
DTP10861	Perio Tip
DTP10862	Universal Tip
DTP10863	Sickle Tip

Blue Air Scaler

Item Number	Description
DTP10870	Sickle Tip
DTP10872	Perio Tip
DTP10875	Universal Tip

P7

Item Number	Description
DTP00690	Beavertail Ultrasonic Scaler Tip
DTP00691	Flat Ultrasonic Scaler Tip
DTP00692	Universal Ultrasonic Scaler Tip
DTP00693	Interproximal Ultrasonic Scaler Tip
DTP00694	Perio Ultrasonic Scaler Tip

LM

Item Number	Description
DTP00784	Universal Tip #39
DTP00785	Flat Tip #32
DTP00786	Special / Perio Tip #33
DTP00788	Endo File Holder Tip 90°
DTP00791	Interproximal #37
DTP00793	Catheter Tip
DTP00794	Scaling tip #31 (Removes Gross Calculus)
DTP00796	Scaling Tip #PE-41PS (Super Slim, Perio)
DTP00791	Set of 5 LM tips, 31, 32, 33, 37, 39



Megasonic

Item Number	Description
DTP100TW	Tip Wrench

Sonix, Lynx, Airstream

Item Number	Description
DTP10882	Universal
DTP10884	Sickle

Sonus V Stack

Item Number	Description
DTP10885	Sonus V Stack with Nose Cone
DTP10886	Universal Tip
DTP10887	Flat Tip

Star Titan

Item Number	Description
DTP10877	Perio Tip
DTP10878	Rebuild Kit

SCALERS & MICROMOTORS

Ram 1000 Micromotor

RAM 1000 Micromotor micro-torque control circuitry design provides smooth constant power at all speeds. It can be operated by hand or foot and goes from a crawl to 35,000 RPM.

- Operates in forward and reverse
- Comes with contra-angle with burs, a diamond disk cutting wheel, and 5 sealed ball bearings (3 in the front nose cone).
- Included Vet Package
- 1 year warranty



Item Number	Description
DTP103RMK	Ram Micromotor

Mystisa Portable Dental Polisher/Micromotor

The brushless motor can perform 10,000 hours or 5 times the life of traditional brush-type motors. It is virtually silent and has no vibration and the sealed motor keeps out all dust and debris. Operates in forward and reverse.

- Digital Color LCD Display Screen
- Battery Level Indicator - Battery Life 18-20 Hours
- Forward / Reverse
- 0-35,000 with Precise Digital Control
- Clear Case with Belt Clip and Handpiece Holder
- Ergonomic Design
- Included Vet Package
- 1 year warranty



Item Number	Description
DTP104RMK	Mystisa Portable Dental Polisher/Micromotor

CONSUMABLES

Prophy Paste

Prophy Cups

Fluoride

Prophy Angles

Contra Heads



PROPHY PASTE

PROPHY ANGLES

CONTRA HEADS

BUR BLOCKS

Prophy Paste

Advanced formula, specifically made for veterinary dentistry. Helps maximize stain removal with minimal enamel loss.

With two different grit options, it makes it easier to tailor the cleaning process to each animal's specific needs when removing biofilm and stains.

- Splatter-free results
- Completely rinses out
- Leaves no residual grit
- Available in Coarse Grit (Medium) and Medium Grit (Fruit)



DTP00730

DTP00729



DTP00731

Item Number	Description
DTP00729	Mint Flavored Propy Paste, Coarse Grit, 9oz
DTP00730	Fruit Flavored Propy Paste, Medium Grit, 9oz
DTP00731	Individual Propy Paste Cups, Fruit Flavor, Package of 100 Medium

Disposable Propy Cups

- Each box contains 144 pieces
- Used to apply propy paste for cleaning and polishing teeth
- Helps effectively remove stains



DTP00720

Item Number	Description
DTP00720	Snap-on Propy Disposable Propy Cups
DTP00670	Screw-on Propy Disposable Propy Cups
DTP00665	Latch-type Propy Disposable Propy Cups

Fluoride

Advanced fluoride formula, in gel or foam form to help protect teeth and discourage bacterial growth. This preventative treatment helps prevent cavities by strengthening the enamel.

- Quick application
- Anti-bacterial effect
- Available in foam or gel



DTP10521

DTP10526

Item Number	Description
DTP10521	Fluoride foam (vanilla), 4.4 oz
DTP10526	Fluoride gel (mint), 4.3 oz

Disclosing Solution

- Used to locate calculus or plaque on teeth
- Great for spot-checking thoroughness of prophylaxis

Item Number	Description
DTP10695	Disclosing Solution 2 oz.



Standard Disposable Propphy Angles

- Individually packed to prevent cross-contamination
- For cleaning and polishing
- High elasticity, effective removal of stains

Item Number	Description
DTP10965	Disposable Propphy Angles (Standard, Pack of 144)



twist™ Propphy Angles

- Individually packed to prevent cross-contamination
- For cleaning and polishing
- High elasticity, effective removal of stains

Item Number	Description
DTP10967	twist™ Propphy Angle (Soft, Bag of 100)
DTP10967-50	twist™ Propphy Angle (Soft, Bag of 50)



Sparkle™ Disposable Prophylaxis Angles

- Beveled gear design runs vibration-free/noise-free
- Flared cup design improves subgingival access and reduces splatter, minimizing cross-contamination
- Light-colored angle with dark colored cup provides greater visibility.
- Single-use

Item Number	Description
DTP10966	Sparkle™ Prophylaxis Angle (Soft, Bag of 100)
DTP10970	Sparkle™ Prophylaxis Angle (Firm, Bag of 100)



DTP10966

True-Seal Prophylaxis Angle

- Sealed head to prevent water and debris from entering the handpiece
- For U-type nose
- Autoclavable and reusable
- 1 year warranty

Item Number	Description
DA10955	True-Seal Prophylaxis Angle, Snap on
DA10956	True-Seal Prophylaxis Angle, Screw-on



DA10955



DA10956

Economy Prophylaxis Angle

- For U-type nose
- Autoclavable and reusable
- 3 month warranty

Item Number	Description
DA551PSN	Economy Prophylaxis Angle, Snap-on
DA552PSC	Economy Prophylaxis Angle, Screw-on



DA552PSC



DA551PSN

Contra Head

- Fits midwest contra-angle attachment
- 14 teeth
- Autoclavable and reusable

Item Number	Description
DAEH-20PK	Snap Type Contra Head
DAEH-20L	Latch Type Contra Head



DAEH-20L



DAEH-20PK

Diamond Burs

- Enhanced cutting performance
- Improved durability
- 5 diamond burs per pack*



Egg

- Preparation of palatal areas
- Occlusal reduction
- Cavity preparation



Item Number	DTP1908C	DTP1908F	DTP390VF
Bur Type	Coarse	Fine	Very Fine
Head Size	018	018	023
Head Length	3.6	3.6	5.3
Total Length	19	19	20

Flame

- Cavity preparation
- Crown and bridge preparation
- *DTP862-012: Pack of 3



Item Number	DTP5128M	DTP5128C	DTP5128F	DTP128VF	DTP862-012*
Bur Type	Medium	Coarse	Fine	Very Fine	Super Coarse
Head Size	012	012	012	012	012
Head Length	8	8	8	8	8
Total Length	21	21	21	21	21

Flat End Cylinder

- Useful for smoothing bone
- Cavity preparation
- Crown and bridge preparation



Item Number	DTP7127C	DTP835-014	DTP837-018	DTP837-014
Bur Type	Coarse	Medium	Coarse Canine	Super Coarse
Head Size	012	014	016	014
Head Length	8	6	8	8
Total Length	21.5	21	21.5	21

Round End Taper


- Crown and bridge preparation
- Round end designed to widen canal wall without the risk of gouging
- *DTP850-016: Pack of 3



Item Number	DTP1610M	DTP1610C	DTP1610F	DTP850-016*
Bur Type	Medium	Coarse	Fine	Super Coarse
Head Size	018	016	016	016
Head Length	10	10	10	10
Total Length	22.5	22.5	22.5	24

Tight Occlusion


- Controlled contouring and finishing
- Smoothing bone
- Trimming and finishing



Item Number	DTP7611TF
Bur Type	12 Flutes FG
Head Size	010
Head Length	4.3

Football


- Preparation of palatinal areas
- Occlusal reduction
- Cavity preparation
- Crown and bridge preparation



Item Number	DTP1923M
Bur Type	Medium
Head Size	012
Head Length	8
Total Length	20


Surgical Length

- Pack of 5
- 25 mm length



Round


- Mark cutting depth in crown preparation
- Used to gain endodontic access
- Prepare cavities



Item Number	DTPFGSL2	DTPFGSL4	DTPFGSL6	DTPFGSL8
Cat. No.	FG. #2	FG. #4	FG. #6	FG. #8
Head Size	010	014	018	023
Head Length	0.8	1.1	1.5	1.9

Cross-Cut Fissure


- Used to section teeth
- Gain endodontic access
- General-purpose bur



Item Number	DTPFGSL557	DTPFGSL558
Cat. No.	FG. #557	FG. #558
Head Size	010	012
Head Length	4.5	4.5

Cross-Cut Tapered Fissure

- Used to section teeth
- Gain endodontic access
- General-purpose bur



Item Number	DTPFGSL700	DTPFGSL701	DTPFGSL702
Cat. No.	FG. #700	FG. #701	FG. #702
Head Size	010	012	016
Head Length	4	4.5	4.5


HP Burs

- Use with low-speed handpieces
- 40 mm
- 5 HP burs per pack



Inverted Cone


- Undercutting cavity preparation
- Produce mechanical interlock



Item Number	DTP10078	DTP10079
Cat. No.	#33 1/2	#35
Head Size	007	010
Head Length	0.7	1.2

Cross-Cut Tapered Fissure


- Used to section teeth
- Gain endodontic access
- General-purpose bur



Item Number	DTP10070	DTP10071	DTP10072
Cat. No.	#700	#701	#702
Head Size	010	012	016
Head Length	4	4.5	4.5

Round

- Mark cutting depth in crown preparation
- Used to gain endodontic access
- Prepare cavities



Item Number	DTP10073	DTP10074	DTP10075	DTP10076	DTP10077
Cat. No.	#1/2	#2	#4	#6	#8
Head Size	007	010	014	018	023
Head Length	0.6	0.8	1.1	1.5	1.9


Latch Type

- Use with low-speed contra-angle latc- type handpieces
- Easily hooks to latch
- 5 latch-type burs per pack



Cross-Cut Tapered Fissure


- Used to section teeth
- Gain endodontic access
- General-purpose bur



Item Number	DTP10036	DTP10038	DTP10042	DTP10039
Cat. No.	#700	#701	#702	#701L
Head Size	010	012	016	012
Head Length	4	4.5	4.5	5.5

Inverted Cone


- Undercutting cavity preparation
- Produce mechanical interlock



Item Number	DTP10040	DTP10067
Cat. No.	#33 1/2	#35
Head Size	007	010
Head Length	0.7	1.2

Round

- Mark cutting depth in crown preparation
- Used to gain endodontic access
- Prepare cavities



Item Number	DTP10059	DTP10060	DTP10064	DTP10066	DTP10068
Cat. No.	#1/2	#2	#4	#6	#8
Head Size	007	010	014	018	023
Head Length	0.6	0.8	1.1	1.5	1.9

Lindemann

- Designed for aggressive cutting of bone structure
- Fine cross cuts and optimum flute depth for efficient cutting and expels waste material
- High-grade fine grain carbide for smooth, vibration-free performance
- Pack of 3



Cross Cut



Item Number	Cat No.	Head Size	Head Length	Total Length
DTP07410	Lindemann162	016	10	25

Plain Cut



Item Number	Cat No.	Head Size	Head Length	Total Length
DTP07411	Lindemann161	016	10	25


Friction Grip (FG)

- 16mm
- 5 FG burs per pack*
- Zekrya 23: Pack of 3, 23 mm length



Cross-Cut Fissure


- Used to section teeth
- Gain endodontic access
- General-purpose bur



Item Number	DTP10021	DTP10019
Cat. No.	#557	#558
Head Size	010	012
Head Length	4.5	4.5

Pear


- Root canal access
- Undercutting dentin



Item Number	DTP10016	DTP10017	DTP10018	DTP10020
Cat. No.	#329	#330	#331	#332
Head Size	006	008	010	012
Head Length	2	2	2	2

Inverted Cone


- Undercutting cavity preparation
- Produce mechanical interlock



Item Number	DTP10033	DTP10034	DTP10035	DTP10037
Cat. No.	#33 1/2	#34	#35	#37
Head Size	007	008	010	014
Head Length	0.7	1	1.2	1.4

Cross-Cut Tapered Fissure

- Section teeth
- Gain endodontic access
- General-purpose bur




Item Number	DTP10025	DTP10026	DTP10028	DTP10030	DTP10032
Cat. No.	#699	#700	#701	#701L	#702
Head Size	009	010	012	016	021
Head Length	4	4	4.5	4.5	5

Item Number	DTP10057
Cat. No.	#699L
Head Size	019
Head Length	5.5

Round

- Mark cutting depth in crown preparation
- Used to gain endodontic access
- Prepare cavities



Item Number	DTP10046	DTP10047	DTP10048	DTP10050	DTP10052
Cat. No.	#1/2	#1	#2	#3	#4
Head Size	007	008	010	012	014
Head Length	0.6	0.6	0.8	1	1.1

Item Number	DTP10054	DTP10056
Cat. No.	#6	#8
Head Size	018	023
Head Length	1.5	1.9

Egg Trimming and Finishing

- 12 flutes: Controlled contouring and finishing



Item Number	DTP07408
Cat. No.	7408
Head Size	023
Head Length	3.5

Zekrya 23 Helical Flute

- Ideal for sectioning an impacted tooth prior to its extraction
- Use for separating roots and removing broken root stump
- Smooth vibration free cutting performance
- Pack of 3*



Item Number	Cat No.	Head Size	Head Length	Total Length
DTP07409	Zekrya23	016	11	23

Perio Scaling

- Scaling and root planing



Item Number	DTP00741	DTP00740
Cat. No.	FG-741	FG-740
Head Size	009	016
Head Length	3	4

Dura-White Stones

- Ideal for finishing and polishing
- Friction grip



Item Number	DTP0450
-------------	---------

Diplomate Bur Kit

Efficient solution for organizing dental burs. This set includes 20 of the most common burs used in veterinary dentistry, all in an autoclavable bur block.

Burs included:

- Round burs (#1/2, #1, #2, #4, #6)
- Cross-cut tapered fissure burs (#700, #701L)
- Two sizes of cylindrical diamond burs (feline, canine)
- Dura-white stones
- 12 extra sections for additional burs




Item Number	Description
DTP10206	Diplomate Bur Block

Magnetic Bur Block with 10 Burs

This set features an autoclavable block and includes a lid (not autoclavable) for convenient storage. It comes with 10 included burs and an additional 4 empty sections to store additional burs.

Burs included:

- Tapered Fissure Cross-Cut (#699, #700, #701, #702)
- Pear (#330, #332)
- Straight Fissure Cross-Cut (#557, #558)
- Round (#6, #8)



Item Number	Description
DTP10201	Magnetic Rectangular Bur Block

Magnetic Bur Block, 7 Bur Capacity

- Autoclavable bur block
- Diameter: 1 - 5/8

Item Number	Description
-------------	-------------

DTP10205	Magnetic Bur Block, 7 Bur Capacity
----------	------------------------------------



Magnetic Bur Block, 28 Bur Capacity

- Autoclavable bur block
- Included lid (not autoclavable)
- Diameter: 2 - 13/14

Item Number	Description
-------------	-------------

DTP10204	Magnetic Bur Block, 28 Bur Capacity
----------	-------------------------------------



Diamond Disc Wheel with Mandrel

- Mounted on 2.35mm shank mandrel
- Disc shape: 358
- Size: 220
- Grit: Coarse
- Thickness: 0.19mm



Item Number	Description
-------------	-------------

DTP10626	Doubled Sided Diamond Disc Wheel with Mandrel
----------	---

Diamond Disc Shield

The disc shield protects operator from spinning diamond cutting disc.

- Designed to fit most straight nose cones

Item Number	Description
-------------	-------------

DTP10627	Diamond Disc Shield
----------	---------------------



HAND INSTRUMENTS

Scalers

Curettes

Excavators

Elevators

Winged Elevators

Periosteal Elevators

Luxators

Levitators

Scissors

Mouth Mirrors

Probes and Explorers

Forceps

Needle Holders

Root Tip Picks

Instrument Trays

Mouth Gags

Retractors

Rodent Instruments

Instrument Kits



WINGED ELEVATORS

PROBES AND EXPLORERS

PETSOFT™

INSTRUMENT KITS

Scalers

These scalers are designed to be used supragingivally (above the gum line).

- Triangular cross-section
- Sharp tip on the end of the blade



Item Number	Description
DTP11315	Diplomat Silk Jacquette Scaler #U15/33
DTP11320	Diplomat Silk Morse Scaler #0-00
DTP11278	Tartar Sickle/Chisel (gold)
DTP11296	Posterior Scaler Petite (red)
DTP11297	Anterior Scaler (red)

DTP11276	Sickle Scaler H6/H7 (green)
DTP11271	Jacquette Feline Scaler 1/2 (black)
DTP11274	Morse Scaler 0/00 (violet)
DTP11275	Sickle Scaler H6/H7 (economy)
DTP11298	Posterior Scaler (red)

HAND INSTRUMENTS | CURETTES

Curettes

- Used to remove subgingival calculus, root planing, and curettage
- One blade with two cutting edges
- Both cutting edges can be used on the front and back of the tooth



DTP11227



DTP11340



DTP11224



DTP11340



DTP11322



DTP00217



DTP00216



DTP11228

Item Number	Description
DTP00216	Diplomate Small Blade Curette-Anterior
DTP00217	Diplomate Small Blade Curette-Posterior
DTP11322	Diplomate Silk Molt Curette #2/3
DTP11228	Barnhart Curette #5/6 (Blue)

DTP11227	Barnhart Curette #5/6
DTP11340	McCalls Curette #17/18
DTP11224	Columbia Curette #13/14

Gracey Curettes

- One cutting edge, designed for deep subgingival calculus, root planing, and curettage of the periodontal pocket.



DTP11226



DTP113112



DTP11313



DTP11225



DTP11695

Item Number	Description
DTP11313	Diplomate Silk Gracey Curette #12/13
DTP113112	Diplomate Silk Gracey Curette #7/8
DTP11695	Gracey Curette #13/14
DTP11225	Gracey Curette #5/6
DTP11226	Gracey Curette #7/8

Diplomate Excavator


Item Number	Description
DTP00211	Diplomate Excavator #31L
DTP00212	Diplomate Excavator #33L Luxator



Regular Tip Elevators

- Blades from 1mm-5mm wide

Item Number	Description
DTP20331	Canine Elevator - Regular
DTP10226	Standard 304 Elevator - Regular
DTP10330	Standard 301 Elevator - Regular



V-Cut Tip Elevator

- Designed to guide it into the socket and cut the periodontal ligament, preventing it from slipping out and causing damage to the gingival tissue

Item Number	Description
DTP10331SM	Feline Elevator w/V-cut - Super Slim




HAND INSTRUMENTS | ELEVATORS

Super Slim Tips

- Excellent for felines

Item Number	Description
DTP10331	Feline Elevator - Super Slim



Serrated Elevators

- Serrations help reduce slippage or rolling off the root of the tooth
- Helps grab the cementum of the sectioned tooth
- Great job of severing the periodontal ligament

Item Number	Description
DTP10338	Serrated Elevator #301
DTP10340	Serrated Mid-Size Elevator #46
DTP20332	Serrated Canine Elevator



DTP10338

Cryer Elevators

- Helpful in removing the caudal teeth in canines

Item Number	Description
DTP10361	Cryer Elevator #30
DTP10362	Cryer Elevator #31



DTP10361



DTP10362

Winged Elevators


Designed by [Dr. Robert Wiggs](#), these are exceptional and aid in secure extractions in both canines and felines.

- Spoon blade allows the instrument to enter the periodontal ligament space
- Sharp tip on the end of the blade
- Generates exceptional leverage between the tooth and surrounding bone when rotated
- Fits well into narrow interdental spaces between teeth


Diplomate Short Shank Winged Elevators

- Semi-circular blade shape adapts to the curved surfaces of roots
- Stainless steel
- Shorter grip

Item Number	Description
DTP10363	Diplomate Winged Elevator 1mm
DTP10364	Diplomate Winged Elevator 2mm
DTP10365	Diplomate Winged Elevator 3mm
DTP10366	Diplomate Winged Elevator 4mm
DTP10367	Diplomate Winged Elevator 5mm
DTP10368	Diplomate Winged Elevator 6mm



DTP10363




DTP10364



DTP10365



DTP10366



DTP10367



DTP10368

HAND INSTRUMENTS | WINGED ELEVATORS

Standard Length

Item Number	Description
DTP10341	Winged Elevators #1 (1mm)
DTP10342	Winged Elevators #2 (2mm)
DTP10343	Winged Elevators #3 (3mm)
DTP10344	Winged Elevators #4 (4mm)
DTP10328	Winged Elevators #6 (6mm)
DTP10329	Winged Elevators #8 (8mm)
DTP10345	Set of 4 Winged Elevators with Case #1-#4



Mini Handle

Item Number	Description
DTP10356	Mini Winged Elevator #1 (1mm)
DTP10357	Mini Winged Elevator #2 (2mm)
DTP10358	Mini Winged Elevator #3 (3mm)
DTP10359	Mini Winged Elevator #4 (4mm)
DTP10360	Set of 4 Mini Winged Elevators with Case #1-#4



Colored with Knurled Handles

- Standard-length elevators with heavy knurled handles
- Knurled to provide a non-slip grip

Item Number	Description
DTP10347	Color Wing Elevator #1 (1mm) Violet
DTP10348	Color Wing Elevator #2 (2mm) Gold
DTP10349	Color Wing Elevator #3 (3mm) Blue
DTP10350	Color Wing Elevator #4 (4mm) Black
DTP10351	Set of 4 Color Wing Elevators with Case #1-#4



Short Shank

Item Number	Description
DTPRX9001S	S/S Short Shank Wing #1
DTPRX9002S	S/S Short Shank Wing #2
DTPRX9003S	S/S Short Shank Wing #3
DTPRX9004S	S/S Short Shank Wing #4
DTPRXSSET	Set of 4 S/S Short Shank #1-4



Standard Length

Item Number	Description
DTPRX9001L	S/S Winged Elevator #1
DTPRX9002L	S/S Winged Elevator #2
DTPRX9003L	S/S Winged Elevator #3
DTPRX9004L	S/S Winged Elevator #4
DTPRXLSET	Set of 4 S/S Winged Elevator #1- 4
DTP00300	Case Only for Winged Elevator



HAND INSTRUMENTS | PERIOSTEAL ELEVATORS

Periosteal Elevators

- Applied once an incision has been made with a scalpel on one or both ends of the diseased pocket to help elevate the gingiva and mucosa away



DTP10335



DTP10339

Item Number	Description
DTP10335	West #2 Perio Elevator
DTP10339	Molt #9 Perio Elevator
DTP10337	Freer #4 Small Perio Elevator



DTP10337

Molt #2/4 Double End

- One end is 2mm and the other is 4mm wide
- Slightly curved shape keeps the instrument on the surface of the bone
- Helps prevent tearing of the raised flap in both canine and feline

Item Number	Description
DTP10333	Molt #2/4 Perio Elevator Double End



Periosteotome Double End

- Most popular periosteal elevator
- Features small, delicate, and decisive straight ends
- Excellent for the gum flaps in both small canines and felines

Item Number	Description
DTP453110	Periosteotome Double End




Luxators

These are specially designed with a fine sharp edge.

- Easy insertion into the sulcular space
- Used to sever the gingival attachment
- Utilized in canines with medium to large-sized teeth
- Reaches segments winged elevators may not

Item Number	Description
DTP0730/1	Luxator Straight 3mm
DTP0730/2	Luxator Straight 5mm
DTP0730/3	Luxator Bent Curve 3mm
DTP0730/4	Luxator Bent Curve 5mm
DTP0730/5	Luxators Set of 4



Diplomate Luxators

Item Number	Description
DTP13101	Diplomate Luxator Curved 1mm, Short Shank
DTP13102	Diplomate Luxator Curved 2mm, Short Shank
DTP13103	Diplomate Luxator Curved 3mm, Short Shank
DTP13104	Diplomate Luxator Curved 5mm, Short Shank



Levitators

The single-ended feline and canine levitators have a "V" cut into the end of the instrument. Easy insertion into the sulcular space

- Helps guide the tip into the socket
- Prevents instrument from slipping and causing damage to gingiva

Curved tips are indispensable for extracting curved-root teeth.

- Designed to follow the concave and convex curvatures of the teeth
- Helps prevent breaking roots

Item Number	Description
DTP00027	Single-ended Canine (Straight)
DTP00028	Single-ended Feline (Straight)
DTP00403	Single-ended Canine (Curved)
DTP00204	Single-ended Feline (Curved)



HAND INSTRUMENTS

Scissors

- Useful in lengthening and trimming the gingival flap

Item Number	Description
DTP20954	Gum and Tissue "S" Shaped Serrated Scissors (4 1/2")
DTP30960	Serrated/Curved Dean Scissors (7")
DTP30961	Diplomate Dean Scissors (6 3/4")
DTP13110	La Grange Scissors (4 3/4")



DTP20954



DTP30960



DTP30961



DTP13110

Mouth Mirrors

- Helps to see around teeth and gum by reflecting an image of the hidden area
- Can be used to direct light onto specific areas

Item Number	Description
DTP11316	Diplomate Silk Mouth Mirror
DTP10545	Dental Mouth Mirror, Size #4



DTP11316



DTP10545

Probes and Explorers

Periodontal probes and explorers are detection instruments.

- Probe is used for measuring the depth of periodontal pockets
- Explorer is used to detect calculus

Item Number	Description
DTP11390	Explorer/Probe #1/23
DTP00012	Perio-probe, Single End #3-6-9-12
DTP00015	Perio-Probe, Single End #3-6-9-12-15
DTP00015	Perio-Probe, Single End #3-6-9-12-15
DTP11291	Probe/Explorer (Blue)



DTP11390



DTP00012



DTP00015



DTP11291

Niemiec Exploration Probe

The color-coded bars assist in concise and visible periodontal pocket diagnostics, with red (6 mm) indicating severity and required advanced treatment.

- Developed by Dr. Niemiec
- Green: If this color is not visible upon a gingival probing during a feline oral exam, a (pathologic) periodontal pocket is present.
- Blue: If this color is not visible upon a gingival probing during a canine oral exam, a (pathologic) periodontal pocket is present.
- Red: Indicates severe problem requiring advanced treatment



Item Number	Description
DTP11324	Niemiec Expro, Color Coded

Extraction Forceps - Right Angle

- Excellent for small incisors and caudal teeth

Item Number	Description
DTP13914	Extraction Forceps - Right Angle (5 1/2")



DTP13914

Fragment Forceps - Angular

- Used for broken root fragments and some small delicate teeth, such as incisors of felines

Item Number	Description
DTP20379	Fragment Forceps - Angular (5")



Tartar Remover Forceps

- Designed to remove heavy deposits of calculus from the teeth

Item Number	Description
DTP11305	Tartar Remover Forceps - Small (4 1/2")
DTP11307	Tartar Remover Forceps - Large (5 1/2")



DTP11305



DTP11307

Rongeur Forceps

- Ideal for small feline teeth
- Break tips are flat and close completely

Item Number	Description
DTP10644	Rongeur Forceps - Small (4 1/2")
DTP20644	Rongeur Forceps - Medium (5 1/2")



DTP10644



DTP20644


Small Forceps

- Excellent for grasping small teeth, especially in felines

Item Number	Description
DTP10645S	Small Extracting Forceps Straight (4 1/2")
DTP10645	Small Extracting Forceps Curved (4 1/2")



DTP10645S



DTP10645

Large Forceps

- Useful for large canine teeth, premolars, and molars

Item Number	Description
DTP20645	Large Extracting Forceps (6 1/2")



DTP20645

Needle Holders

Item Number	Description
DTP11325	Castroviejo Needle Holder w/ Scissor
DTP11327	Olsen-Hegar STR T/C 15cm



DTP11325

Root Tip Picks

- Useful for removing small pieces of broken root tips

Item Number	Description
DTP10175	Root Tip Pick #78 - Straight
DTP10176	Root Tip Pick #79 - Left Angle
DTP10177	Root Tip Pick #80 - Right Angle
DTP10178	Root Tip Pick - Set of 3



DTP10175



DTP10176



DTP10177

Instrument Trays

- Enhance instrument organization

Item Number	Description
DTP10925L	Large Instrument Tray (Flat) 13 1/2" x 9 5/8"
DTP10925D	Large Instrument Tray (Divided) 13 1/2" x 9 5/8"
A01-115	Stainless Steel Tray 13 1/2" x 9 5/8"
DTP10925	Small Instrument Tray (Divided) 8" x 4"



DTP10925D



A01-115



DTP10925



DTP10925L

Mouth Gag

- Do not autoclave

Item Number	Description
DTP42254	Large Mouth Gag, 7"
DTP42250	Canine Mouth Gag, 5"
DTP42256	Feline Mouth Gag, 3"



DTP42254



DTP42250



DTP42256

Plastic Cheek Retractors

- Ideal for pulling back the lip corners of canines or felines for access to the caudal part of the mouth

Item Number	Description
DTP42257	Small, 1 1/2" wide - feline
DTP32590	Large, 2" wide - canine



DTP42257



DTP32590

Senn Retractor

- Indispensable for gingiva and cheek retraction
- Right-angled, single blade at one end
- Curled, three-pronged claw at the other end

Item Number	Description
DTP42260	Senn Retractor



Atraumatic Lip Retractor


- Developed by Dr. Tony Woodward
- Gentle and flexible
- Ends of the instrument include jaws that can be gently bent to accommodate different lip thicknesses
- Do not autoclave

Item Number	Description
DTP42258	Atraumatic Lip Retractor




PetSoft™ Home

- Single multi-use for home care designed with pet owners and clinics in mind
- Safe and effective retractor used to look and monitor dental health or injuries that occur in the mouth at home
- Ergonomic, soft, and bendable tip offering comfort to companion animal
- Easily monitor tartar build-up and view the entire mouth
- Simple to clean and keep hygienic by rinsing with warm water after use
- 12 large and 12 small sizes per pack
- Dishwasher safe
- Dimensions (Small): Weight: 0.40 oz; Overall Size: 18 cm long, 1.4 cm handle size; Delta Shape Size: 2.6 cm long, 2.7 cm wide, 0.53cm thick
- Dimensions (Large): Weight: 0.50oz; Overall Size: 18.7 cm long, 1.4 cm handle size; Delta Shape Size: 3 cm long, 3.1 cm wide, 0.53 cm thick





Item Number	Description
DTP10973	PetSoft Home, 24 Pack for Home Care



PetSoft™ Pro

- Designed exclusively for professional use at clinics to improve access, visibility, and comfort during procedures and exams
- Ergonomic, soft, and bendable tip to customize retraction
- Textured surface to prevent slippage of tongue and other soft tissues
- Non-reflective material for laser procedures
- Autoclavable for 30 cycles
- 6 large and 6 small sizes per pack
- Dimensions (Small): Weight: 0.40 oz; Overall Size: 18 cm long, 1.4 cm handle size; Delta Shape Size: 2.6 cm long, 2.7 cm wide, 0.53cm thick
- Dimensions (Large): Weight: 0.50oz; Overall Size: 18.7 cm long, 1.4 cm handle size; Delta Shape Size: 3 cm long, 3.1 cm wide, 0.53 cm thick


Item Number	Description
DTP10972	PetSoft Pro (Autoclavable), 12 Pack for Professional



Molar Extraction Forceps

- Designed to grip without crushing
- Most abscesses around the rabbit head have their origin in molar tooth roots

Item Number	Description
DTP454600	Molar Extraction Forceps



Molar Rasp

- Used to remove tooth spurs without damaging the delicate root structure


Item Number	Description
DTP454800	Molar Rasp



Rodent Molar Cutters

- Used to trim spurs and hook angles
- Unique shape enables the surgeon to see what is being cut

Item Number	Description
DTP454650	Rodent Molar Cutters



Rodent Bur Guard

- Soft tissue protector used on low-speed handpieces

Item Number	Description
DTP02443	Rodent Bur Guard



Rabbit Incisor Luxator

Item Number	Description
DTP454700	Rabbit Incisor Luxator



Rodent Mouth Gag

- Essential tool when examining the mouth of a rabbit

Item Number	Description
DTP454500	Rodent Mouth Gag




Cheek Dilators

- Allow visualization of molars and protects cheeks from instruments
- Useful as a gag for guinea pigs and chinchillas

Item Number	Description
DTP454550	Small Rodent Cheek Dilator
DTP454555	Large Rodent Cheek Dilator



DTP454550



DTP454555

Rodent Dental Instrument Kit

- Designed specifically for rodent patients

Item Number	Description
DTP454850	Rodent Dental Kit

This set includes the following instruments:


DTP00300	Instrument Case
DTP454500	Rodent Mouth Gag
DTP454550	Small Rodent Cheek Dilator
DTP454555	Large Rodent Cheek Dilator
DTP454600	Rodent Molar Extractor
DTP454650	Rodent Molar Cutter
DTP454700	Rabbit Incisor Luxator
DTP454800	Molar Rasp



DTP454500




DTP454550



DTP454555



DTP454800



DTP454650



DTP454600



DTP454700

Diplomate Perio Surgery Kit

- Autoclavable stainless steel cassette

Item Number	Description
-------------	-------------

DTP4-008122	Diplomate Perio Surgery Kit
-------------	-----------------------------

This set includes the following instruments:

DTP00211	Excavator #31L
DTP00212	Excavator #33L
DTP00216	Small Blade Curette-anterior
DTP00217	Small Blade Curette-posterior
DTP11313	Silk Gracey Curette #12/13 1
DTP113112	Silk Gracey Curette #7/8
DTP11314	Silk Expro 23E/18P
DTP11315	Silk Jacquette Scaler #U15/33
DTP11316	Silk Mouth Mirror
DTP11320	Silk Morse Scaler #0-00
DTP11322	Silk Molt Curette #2/3
DTP11325	Needle Holder
DTP13107	Adson-Brown Tissue Forceps
DTP13108	Scalpel Blade #3
DTP13110	La Grange Scissors
DTP453110	Periosteotome
DTP4008122	Autoclavable Cassette



Diplomate Extraction Kit

- Autoclavable stainless steel cassette

Item Number	Description
-------------	-------------

DTP3-188114	Diplomate Extraction Kit
-------------	--------------------------

This set includes the following instruments:

DTP10179	Woodward Heidbrink Root Tip
DTP10363	Diplomate Winged Elevator 1 mm
DTP10364	Diplomate Winged Elevator 2 mm
DTP10365	Diplomate Winged Elevator 3 mm
DTP10366	Diplomate Winged Elevator 4 mm
DTP10367	Diplomate Winged Elevator 5 mm
DTP10368	Diplomate Winged Elevator 6 mm
DTP10646	Diplomate Friedman Rongeur 4.5
DTP11314	Silk Expro 23e/18p
DTP11322	Silk Molt Curette #2/3
DTP13103	Luxator Curved 3 mm
DTP13104	Luxator Curved 5 mm
DTP13107	Adson-Brown Tissue Forceps
DTP13108	Scalpel Blade #3
DTP13109	Olsen-Hegar Needle Holder
DTP30961	Diplomate Dean Scissors
DTP00212	Excavator #33L
DTP4008122	Autoclavable Cassette




Diplomate Extraction Pack

- Fundamental dental instruments required for extractions in canines and felines

Item Number	Description
DTP10616	Diplomate Diplomate Extraction Pack

This set includes the following instruments:

DTP13101	Luxator Curved 1 mm
DTP13102	Luxator Curved 2 mm
DTP13103	Luxator Curved 3 mm
DTP13104	Luxator Curved 5 mm
DTP13105	Small Breed Extracting Forcep
DTP13106	Baby Molt Periosteal Elevator
DTP13107	Adson-Brown Tissue Forcep
DTP13108	Scalpel Blade #3
DTP13109	Olsen-Hegar Needle Holder
DTP11313	Silk Gracey Currette #12/13
DTP4008122	Cassette for Diplomate Extraction Kit




Diplomate Perio Pack

- Ideal for periodontal scaling, calculus debridement, and root planing

Item Number	Description
DTP10618	DentalAire Perio Pack

This set includes the following instruments:

DTP11311	Silk Gracey Currette #1/2
DTP11312	Silk Gracey Currette #7/8
DTP11313	Silk Gracey Currette #12/13
DTP11314	Silk Expro 23E/18P
DTP11315	Silk Jacquette Scaler #U15/33
DTP00300	Case
DTP11316	Silk Mouth Mirror



HAND INSTRUMENTS | KITS

Woodward Extraction Kit

Designed by Dr. Tony Woodward, the Woodward Extraction Kit includes all fundamental dental instruments required for extractions in canines and felines.

Item Number	Description
DTP10613	Woodward Extraction Kit

This set includes the following instruments:

DTP10341	S/S Winged Elevator 1mm
DTP10342	S/S Winged Elevator 2mm
DTP10343	S/S Winged Elevator 3mm
DTP10344	S/S Winged Elevator 4mm
DTP0730/3	Curved Luxator 3mm
DTP0730/4	Curved Luxator 5mm
DTP30960	Curved Serrated Dean Scissors
DTP10333	Molt 2-4 Perio Elevator
DTP10644	Small Rongeur Forceps, 4.5"
DTP10175	Straight Root Tip Pick



Mini Short Shank Hand Instrument Kit

- Ideal for feline dentistry and scaling teeth of rodents

Item Number	Description
DTP11287	Mini Short Shank Hand Instrument Kit

This set includes the following instruments:

DTP11282	M5/6 Gracey #5/6
DTP11283	M2/3 Jacquette #2/3
DTP11284	M0/00 Morse #0/00
DTP11285	M6/7 Sickle Scaler H6/H7
DTP11286	M5/6 Barnhart #5/6
DTP11288	M15/16 Tartar Sickle Chisel
DTP10925	Instrument Tray



Basic Dental Instrument Set-up Kit

Item Number	Description
-------------	-------------

DTP10612	Basic Dental Instrument Set-Up Kit
----------	------------------------------------

This set includes the following instruments:

DTP10545	Mouth Mirror
DTP11227	5/6 Barnhart Curette Eco
DTP11275	H6/H7 Scaler Eco
DTP11279	Tartar Remover Scraper
DTP11390	Probe / Explorer 1/23
DTP11226	Gracey Curette 7/8 Eco
DTP10925	Instrument Tray



Color Coded Instrument Kit

- Set of 6 instruments
- Autoclavable tray to easily organize and clean instruments

Item Number	Description
-------------	-------------

DTP11330	Set of 6 with Autoclavable Tray
----------	---------------------------------


This set includes the following instruments:

DTP11225	Gracey 5/6 Curette
DTP11228	Barnhart 5/6 Curette C/C
DTP11271	Feline Scaler
DTP11274	Morse Scaler C/C
DTP11276	H6/H7 Scaler
DTP11278	Tartar Scraper / Chisel C/C
DTP10925	Small Animal Instrument Tray



SimpliScale Periodontal Kit


- Exclusive universal/area-specific designed instrument set
- Addresses all periodontal needs for clinicians
- Unique instruments allow clinicians to eliminate the need to change tools while performing lingual and buccal cleaning techniques
- Ergonomic non-slip grip
- Wear resistant DuraGradeMax steel with a lightweight alloy core
- Autoclavable instruments and cassette
- Non-fade color-coded instrument handles



Item Number	Description
DTP11295	SimpliScale Periodontal Kit

This set includes the following instruments:

DTP11294	LM DentalAire Universal Gracey Regular
DTP11292	Universal Gracey Petite
DTP11293	Universal Gracey Anterior
DTP11298	LM DentalAire Scaler Posterior
DTP11296	Scaler Posterior Petite
DTP11297	LM DentalAire Scaler Anterior
DTP11291	LM DentalAire Probe/Explorer
DTP11299	LM 8 Instrument Cassette
DTP90605	Acrylic Test Stick Pack of 3



ADVANCED DENTAL CARE

Bonded Sealant Kit

Curing Light

Restorative

Pluggers

K-reamers

H-files

Gutta Percha Points

Absorbant Paper Points

Sharpening Instruments

Impression Products



BONDED SEALANT KIT

DENTALAIRE

ITEM # DTP99950



CURING LIGHT

PLUGGER/SPREADERS

GUTTA PERCHA POINTS

SHARPENING KIT

Bonded Sealant Kit

The **Bonded Sealant Kit** is used to help smooth the fracture site, decrease plaque adherence, calculus accumulation, and decrease soft tissue irritation from the fractured surface. Plugging the dentin tubules with the bonded sealants minimizes sensitivity and helps prevent infection of the pulp chamber.

This kit has all you need to start applying bonded sealants for exposed vital dentin and bonding composite restorations and can treat about 30-44 teeth.

Kit Includes:

- Etch 37% (5mg Syringe) - Color Blue
- One-Step Universal Light Cure Dental Adhesive
- Fortify Surface Sealant
- Super-Snap Mini Kit
- Disposable Mixing Wells (48 count)
- Disposable Applicators (Orange, 100 count)
- Blue Tip Etch Gel Pack (10 count)
- Latch Type Contra Angle Eco.



Item Number	Description
DTP99950	Bonded Sealant Kit
DTP99955	Restorative Bonded Sealant Kit #1 with Clear Blue Curing Light
DTP99956	Restorative Bonded Sealant Kit #2 with Clear Blue Curing Light, Flowable Composite, general-purpose composite, and the Woodson PFI Plastic Filling Instrument

Clear Blue Curing Light

The Clear Blue LED Curing Light is portable, quiet, and powerful. The digital display on the handle is simple to read and offers 3 different modes (Strong, Soft Up, and Pulse) for the device.

- With LED light, it never heats up and offers a power supply of 100V to 240V
- Curing times are faster than traditional halogen lights
- Kit includes: handle, light guide rod, charger, and charger base
- 1 year warranty excludes parts and accessories

Item Number	Description
DTP00683	Clear Blue Curing Light
DTP00689	Angled Light Guide
DTP00683-C	Clear Blue Curing Light Charging Base
DTP00683-B	Clear Blue Curing Light Battery



DTP00689



DTP00683

Super-Snap® Polishing System

Designed for easy, safe, and fast contouring, finishing, and polishing of both micro-filled and hybrid composites.

Kit Includes:

- 180 disks (100 standard, 80 mini)
- 40 poly strips
- 2 dura-white stones (shapes CN1 and FL2)
- 1 Composite Fine and Mid-Point
- 4 CA Mandrels



Item Number	Description
DTP24490	Super-Snap® Polishing System

Super-Snap® Mini Kit

Kit Includes:

- 48 disks (32 standard, 16 mini)
- 1 Composite Fine and Mid-Point
- 1 CA mandrel



Item Number	Description
DTP24491	Super-Snap® Mini Kit

Replacement Disks

Item Number	Description		
DTPS230747	Black Mini Replacement Disks (Coarse, Pack of 50)	DTPS230815	Red Mini Replacement Disks (Super Fine, Pack of 50)
DTPS230805	Green Mini Replacement Disks (Fine, Pack of 50)	DTPS230820	Violet Mini Replacement Disks (Medium, Pack of 50)

Dura-White Stones

- Ideal for finishing and polishing
- Friction grip

Item Number	Description
DTP0450	Dura-white Stone FG



Vitrebond Glass Ionomer Liner/Base

- Dual-cured (both light and self-cure) glass ionomer
- Useful as a base liner in endodontic restoration and restoration of caries at gingival margin
- Fluoride release aids in preventing the extension of carious lesions

Item Number	Description
DTP24492	Vitrebond Glass Ionomer Liner/Base



One-Step®

- Light-cured single-component adhesive
- Bonds to a multitude of dental substrates, including self-cured and dual-cured materials

Item Number	Description
DTP1102P	One-step 6ml bottle



One-Step Plus®

- Light-cured single-component adhesive
- Bonds to a multitude of dental substrates, including self-cured and dual-cured materials
- Formulated with filler particles to help reduce post-operative sensitivity



Item Number	Description
DTP2001P	One-step Plus 6ml Bottle

All-Bond 2® Primer A & B Set

- Dual-cured adhesive material designed for use in all bonding procedures
- Comprehensive ability to bond dentin, enamel, composite, precious, semi-precious, and non-precious casting alloys, and silane-treated porcelain



Item Number	Description
DTP2510A	All-bond 2 Primer A & B Set

All-Bond 2® Kit with Uni-Etch® 32% w/BAC

Universal dental adhesive, which is still unsurpassed in its ability to bond to dentin, enamel, new or old composites, and amalgams. It has a dual-cured primer, that can be used with or without a light cure unit. Includes an antimicrobial etchant that disinfects at the same time.

Kit Includes:

- Primers A & B (6ml bottle)
- D/E Bonding Resin (6ml bottle)
- Pre-Bond Resin (3ml bottle)
- Syringe
- Syringe Uni-Etch w/BAC (5g)



Item Number	Description
DTP42500	All-Bond 2 w/Etch, DE Resin, Prebond Kit

All-Bond 3® Bottle Kit

- Micro-mechanically bonds to substrates
- Compatible with light, dual, and self-cured materials
- Offers increased hydrophobicity for a long-lasting adhesive bond
- HEMA-free, radiopaque
- Reduces chances of misdiagnosing caries



Item Number	Description
DTP36200K	All-bond 3 Bottle Kit

All-Bond Universal® Bottle

- Not moisture sensitive
- Use on wet, dry, or moist tooth surface
- Virtually no post-operative sensitivity
- Contains MDP monomers for enhanced durability
- Bonds to all substrates



Item Number	Description
DTP7202P	All-bond Universal 6ml

Uni-Etch® 32% Syringe

- 32% semi-gel phosphoric acid etchant
- Blue color for easy visualization
- Etchant rinses away cleanly and quickly
- Creates micro retentive surfaces that are necessary for successful bonding



Item Number	Description
DTP5501E	Uni-Etch 32% 5gm Syringe

Etch-37™ w/BAC

- 37% semi-gel phosphoric acid etchant
- Blue color for easy visualization
- Etchant rinses away cleanly and quickly
- Creates micro retentive surfaces that are necessary for successful bonding

Item Number	Description
DTP5503E	Etch-37™ 37% 5gm Syringe
DTP5638E	Etch Gel 37% 30gm Bulk Bottle



Phosphoric Etchant Syringe Tips

- Light Blue Phosphoric Etchant Syringe Tips used with Uni-Etch® and Etch-37™
- 30 per pack

Item Number	Description
DTP5503F	Syringe Tip



Fortify™

- Light-cured, low-viscosity resin formulated for use as a composite surface sealant
- Improved wear resistance and long-term marginal integrity
- Effective with all anterior/posterior composites

Item Number	Description
DTP9031F	Fortify 5ml bottle



Mixing Wells

- Ideal for adhesives, opaques and all materials dispensed in small quantities

Item Number	Description
DTP80496	Mixing Wells



Applicators

- Applies a greater and more uniform amount of material
- Disposable applicators for applying adhesives
- Pack of 100

Item Number	Description
DTP80247	Blue Bristle Brush Applicator
DTP80241	Regular Orange Applicators



Mixing Pads

- Great to use with Vitrebond
- 5 per pack
- 2 3/8" x 2 3/8"

Item Number	Description
DTP80485	Mixing Pads



Flour Pumice

- Fine grit
- 8 oz

Item Number	Description
DTP10860	Pumice 8oz



Aeliteflo™ Syringe Shades

These flowable composites are versatile and can be dispensed and injected into the preparation without dripping or slumping.

- Easy handling for effortless adaption to dentition
- Micro-hybrid formula provides outstanding strength
- Ideal for Class V restorations, Class III's and small Class IV's

Item Number	Description
DTP671A2	Shade A-2, refill 1.5gm. Syringe, Light Cure/Flowable
DTP671A3	Shade A-3, refill 1.5gm. Syringe, Light Cure/Flowable
DTP673C2	Shade C-2, refill 1.5gm. Syringe, Light Cure/Flowable
DTP674D3	Shade D-3, refill 1.5gm. Syringe, Light Cure/Flowable



DTP671A2

Aelite™ All-Purpose Body Syringe

Aelite All-Purpose Body is a durable, light-cured, universal micro-hybrid composite suitable for all applications from Class I through Class VI.

- High compressive strength
- Designed to blend with surrounding dentition
- Radiopaque for easy detection on radiographs
- 4g syringe

Item Number	Description
DTP814A2	Aelite All-Purpose Body Composite (Shade A-2)
DTP814A3	Aelite All-Purpose Body Composite (Shade A-3)
DTP814C2	Aelite All-Purpose Body Composite (Shade C-2)
DTP814D3	Aelite All-Purpose Body Composite (Shade D-3)



DTP814A2

Visible Light Cure Calcium Hydroxide Liner

- Used as a liner in cases of near pulp exposure and for vital pulp therapy after placement of calcium hydroxide powder on the exposed pulp tissues
- Light cure polymerization is convenient and shortens clinical procedure times

Item Number	Description
DTP10575	Visible Light Cure Calcium Hydroxide Liner

Protemp Plus Intro Kit

Protemp Plus is a composite material for the fabrication of temporary restorations. Based on multi-functional methacrylate esters, this two-component system has filling material-like properties.

Kit Includes:

- Base paste/Catalyst Paste (1 count, 67g)
- Garant Dispenser 10:1 (1 count)
- Garant Mixing Tips - Blue (16 count)
- RelyX Temp NE Base Paste (1 count, 30g)
- Catalyst Paste (1 count, 13g)
- Filtek Supreme Ultra Flowable Paste A 2 (2 count, 4g)



Item Number	Description
DTP46201	Protemp Plus Intro Kit

Woodson PFI #3 Plastic Filling Instrument

- Useful for placing and contouring composite restorative materials

Item Number	Description
DTP00210	Woodson PFI #3 Plastic Filling Instrument



Cement Spatula

Item Number	Description
DTP10212	Cement Spatula



Zinc Oxide Powder

- Formulates cement and root canal sealers

Item Number	Description
DTP10230	Zinc Oxide 1.5oz, Powder



Holstrom Plugger/Spreaders

- Designed by Dr. Steve Holmstrom
- Useful for obstructions in larger canine teeth

Item Number	Description
DTP453120	Holmstrom Plugger/Spreader #20
DTP453135	Holmstrom Plugger/Spreader #35
DTP453150	Holmstrom Plugger/Spreader #50
DTP453165	Holmstrom Plugger/Spreader #65
DTP453190	Holmstrom Plugger/Spreader #90
DTP453191	Holmstrom Plugger/Spreader Set of 5 (#20, #35, #50, #65, #90)



Diplomate Pluggers

- Designed for multiple applications, including compaction of gutta-percha (when heated), placement of composite restoratives, treatment of medicaments when performing vital pulp therapy (pulp caps), and apical endodontic restorations

Item Number	Description
DTP00016	Diplomate Plugger A, .03mm
DTP00017	Diplomate Plugger B, .04mm
DTP00018	Diplomate Plugger C, .05mm
DTP00019	Diplomate Plugger D, .06mm
DTP00020	Woodward Pluggar Set (A, B, C, D)



Finger Pluggers, 25mm

- Used to vertically condense gutta percha in root canal procedures
- Set of 6 assorted, 25mm, color-coded
- Sizes #15 - 40

Item Number	Description
DTP190662	Finger Pluggers 25mm



K-reamers, 60mm

- Set of 3 assorted
- 60mm

Item Number	Description
-------------	-------------

DTP01777	#15-20-25
DTP01778	#30-35-40
DTP01779	#45-50-55
DTP01780	#60-70-80



H-Files, 120mm

- Set of 3 assorted
- 120mm

Item Number	Description
-------------	-------------

DTP01761	#20-25-30
DTP01762	#35-40-45
DTP01763	#60-70-80



H-Files, 60mm

- Set of 3 assorted
- 60mm

Item Number	Description
-------------	-------------

DTP01747	#15-20-25
DTP01748	#30-35-40
DTP01749	#45-50-55
DTP01750	#60-70-80




Color-coded Barbed Broach

- Used to remove root canal contents
- 21mm
- Set of 6

Item Number	Description
-------------	-------------

DTP711330	Barbed Broach, 21mm
-----------	---------------------



Eugenol Liquid

Non-irritant liquid to be mixed with Zinc Oxide powder in order to fabricate a low-strength Zinc Oxide Eugenol (ZOE) temporary cement.

- Serves as a temporary tooth filling or as a base cement to affix a temporary tooth filling
- Affix dental devices such as crowns or bridges
- Applied to a tooth to protect the tooth pulp



Item Number	Description
DTP10153	Eugenol Liquid 1 oz.

Finger Spreaders

- Used for laterally spreading gutta percha in root canal procedures
- Set of 6 assorted, 25mm, color-coded
- Sizes #15-40



Item Number	Description
DTP190676	Finger Spreaders 25mm

Gutta-Percha Points

- Designed for animal canals
- Oversized
- 60mm
- 60 pieces




Item Number	Description	Item Number	Description
DTP122205	#15, #20, #25 (20 of each size)	DTP122219	#130
DTP122206	#30, #35, #40 (20 of each size)	DTP122220	#140
DTP122211	#45, #50, #55 (20 of each size)		
DTP122214	#60, #70, #80 (20 of each size)		
DTP122215	#90		
DTP122216	#100		
DTP122217	#110		
DTP122218	#120		



Absorbant Paper Points

- Designed for animal canals
- Oversized
- 60mm
- 60 pieces



Item Number	Description
DTP218291	#15, #20, #25 (20 of each size)
DTP218292	#30, #35, #40 (20 of each size)
DTP218293	#45, #50, #55 (20 of each size)
DTP218294	#60, #70, #80 (20 of each size)
DTP218295	#90
DTP218296	#100
DTP218297	#110
DTP218298	#120
DTP218299	#130
DTP218300	#140

Endo Locking Pliers

- Useful for handling paper points and gutta-percha points

Item Number	Description
DTP10570	Endo Locking Pliers



Sharpening Instruments

Item Number	Description
DTP10615	Arkansas Sharpening Stone, Conical 3" base x 5/16"
DTP10375	Arkansas Sharpening Stone, Flat 4" x 1" x 1/4"
DTP10660	Lubricating Oil, 0.5 oz.
DTP90605	Acrylic Test Sticks (Pack of 2)
DTP10616K	Sharpening Stone Kit (includes items above)



DTP10615



DTP10375



DTP10660



DTP90605

Alginate Powder

Item Number	Description
DTP10125	Alginate Powder



Alginate Spatula

Item Number	Description
DTP10213	Alginate Spatula



Impression Tray

- Set of 8

Item Number	Description
DTP00966	Impression Tray



Flexi-Bowl

- 700ml
- Color: dark green

Item Number	Description
DTP150110	Flexi-Bowl



Protective Gear

- Splash Shield: protects from paste splatter and bacteria; included band

Item Number	Description
DTP11234	Splash Shield, Pack of 3
DTP360005	Safety Glasses



DTP11234



DTP360005

EDUCATIONAL RESOURCES

Mouth Models

Dental Charts

Fig. 5

Fig. 6

Non-vital teeth should always be treated. Many discolored teeth are non-vital (dead). The beige color is non-vital. A dental radiograph of this tooth (fig. 6) shows that the root canal system is much wider than the

Fig. 8

Fig. 9

Some non-vital teeth are not treated. Many primary teeth are non-vital (fig. 9) and may be treated with a pulpotomy (fig. 8). The wider root canal system commonly reveals pathology in teeth that



MOUTH MODELS


SKULLS

DENTAL CHARTS

DENTAL POSTER

Mouth Model

Item Number	Description
DTP101VMD	Canine
DTP102VMC	Feline
DTP103VMDC	Combo Canine and Feline



DTP102VMC

DTP101VMD

Mouth Model with X-ray Teeth

- Radio-opaque natural removable and replaceable teeth
- Transparent, firm, and smooth

Item Number	Description
DTP105VMC	Feline Mouth Model w/X-ray Teeth
DTP105VMD	Canine Mouth Model w/X-ray Teeth



DTP105VMC

DTP105VMD

Mouth Model with Gingiva

Item Number	Description
DTP101VMDCO	Canine with gingiva
DTP102VMDCCO	Feline with gingiva



DTP102VMDCCO

DTP101VMDCO

Skulls

Item Number	Description
DTP62200	Canine Skull
DTP62300	Feline Skull




DTP62300

DTP62200

Stick on Veterinary Dental Charts

Item Number	Description
DTP20340	Starter Kit - 25 Canine / 25 Feline
DTP20341	Canine Refill - 50
DTP20342	Feline Refill - 50



Dental Pathology Poster

Item Number	Description
DTP61002	DentalAire Dental Pathology Poster "Dental Problems You See Everyday"



Radiology DVD

Item Number	Description
DTP62500	Radiology Simplified Instruction Educational DVD by Dr. Tony Woodward and Dr. Brook Niemiec Includes Dental Radiology Positioning Instructions

LIMITED WARRANTY

Dentalaire warrants its products to be free from defects in materials and workmanship only. No other warranties are expressed or implied.

Air Compressors and Portable Units

All compressors and portable dental units including wall-mounted and countertop units are warranted for one (1) year with the exception of:

- Prestige Dental Station (2 years)
- Diplomat Dental Station (2 years)

X-ray Imaging

- DTX Digital Dental Imaging Sensors (3 years)
- DC X-ray Generator Wall Mount and Mobile (2 years)
- DHX Handheld X-ray Generator (1 year)
- DTX IO Camera (1 year)
- DHX Handheld Battery (6 months)

Miscellaneous Products

- Water Distiller (1 year)
- Service Parts (3 months)

Instruments

- Dentalaire and VetCare instruments (1 year)

If any of our instruments break, chip, or bend during any dental procedure, Dentalaire will replace the instruments at no charge. Defective item(s) must be returned to us for exchange.

Terms and Conditions

Terms of payment will be established by results of credit report/references. Once credit terms are established, invoices will reflect terms extended by Dentalaire:

- | | |
|--|-------------------------|
| • Order placement/information to include | • Product Description |
| • Invoice to/Ship to information | • Quantity |
| • Purchase Order Number | • Desired Shipping Date |
| • Product Number | • Method of Shipment |

Note: Opened handpieces and equipment may not be returned for credit, but will be repaired or replaced in accordance with manufacturer warranties. All warranty items, anything after 90 days, customer is responsible for the freight. Please inspect all packages carefully prior to opening.

Returns

All items must be returned within 30 days of purchase in like new condition with all original packaging. A 25% restocking fee will apply to all returns. Returns will be inspected before a credit is issued. Items cannot be returned after 30 days.

All refurbished items are final sale and not returnable.

Disclaimer

THE ABOVE WARRANTIES ARE VOID IF ITEMS ARE CARELESSLY USED, IMPROPERLY INSTALLED, IMPROPERLY STERILIZED, AND/OR MISUSED. NO CLAIM FOR LABOR OR CONSEQUENTIAL DAMAGES WILL BE ALLOWED. SHIPPING IS THE RESPONSIBILITY OF THE CUSTOMER. ANY DISCREPANCIES IN YOUR ORDER MUST BE REPORTED WITHIN 10 DAYS OF RECEIPT OF THE ORDER.

Handpieces and Accessories on Dental Stations

Proof of invoice is required and defective scaler tips are to be sent back:

- High-speed handpiece (1 year)
- High-speed turbine (3 months)
- Low-speed handpiece (1 year)
- True Seal prophyl angles (1 year)
- Economy prophyl angles (3 months)
- Scaler tips (3 months)
- LM Piezo Ultrasonic Scaler handpiece (1 year)
- Metal Piezo Handpiece (1 year)

All other prophyl angles and tips are not applied. Warranty is limited repair, parts, or replacement of item by Dentalaire. All other warranties are described in the individual product description.

Scalers and Micromotors

- Newtron® Booster (2 years, 90 days scaler tips)
- Newtron® P5 with LED (2 years, 90 days scaler tips)
- Piezotome® Cube (2 years, 90 days scaler tips)
- Cocoon Spray (1 year, 90 days scaler tips)




Parts

Need parts for your Dentalaire equipment? We offer a wide selection of parts to restore, repair, and upgrade your products.

Visit our website to find the part you need, or contact our technical support specialist for expert guidance.

www.dentalaireproducts.com/parts



Unlimited Customer Service and Technical Support

Dentalaire offers unlimited customer support. From troubleshooting dental stations and scheduling software demos to repairs, our representatives are ready to assist you.

Contact us today to learn more.

Call: (800)866-6881

E-mail: customerservice@dentalaire.com

Website: www.dentalaireproducts.com/contact

Business Hours: 7:30am - 4pm PST (Monday - Friday)



Scan to Visit Our Website



17742 Mitchell North STE B, Irvine, CA, 92614 | info@dentalaire.com

(800) 866-6881 | www.dentalaire.com