




DENTALAIRE

Meeting and Exceeding the Needs of Veterinary Dentistry since 1988

**Portable
X-Ray**

Wireless Portable X-Ray

Anytime, Anyplace.



**Handy and Powerful
Excellent Image
Ergonomic Design**

www.Dentalaire.com | 800.866.6881

DentalAire | 17742 Mitchell North, Suite B Irvine CA 92614

info@dentalaire.com



DENTALAIRE


Meeting and Exceeding the Needs of Veterinary Dentistry since 1988

Portable X-Ray



Total weight 2Kg, 3mA X-ray Tube current

- Simple operation
- Wide graphic LCD
- Lightweight and ergonomic design
 - Highest image quality and lowest patient dosage
 - Zero radiation leakage
 - High frequency, constant-potential X-ray generator
- Optimized to both analog film and digital sensor



Technical Spec's

X-Ray Generator	DC High Frequency	Biox
X-Ray Insulation	Tube Voltage	60kV
X-Ray Tube	Tube Current	3mA
	Anode Type	Stationary
	Anode Angle	20°
	Focal spot size	0.8mm
	Heat capacity	8.5Khu
	0.01~1.6 sec(0.01 step)	
Exposure Time	1.6mm Al	
Total Filtration	100mm, 210mm(long cone)	
SSD	2.2kg	
Total Weight	180(w) x 73(d) x 138(h) mm	
Dimension	Battery	DC 24V
Power Requirement	Battery Charger	AC 100~240V
	Power Supply(Stand)	AC 100~240V

www.Dentalaire.com | 800.866.6881

Dentalaire | 17742 Mitchell North, Suite B Irvine CA 92614

info@dentalaire.com