


Countertop III Dental Station



DTP00504

The DentalAire Countertop Unit is compact, powerful, and will easily integrate into your veterinary dentistry practice.

This countertop unit includes:

- High-speed Push Button Drill or LED High-speed Handpiece
- Low-speed Polisher Handpiece
- LM Piezo Scaler with one #39 Tip Universal and #32 Flat Tip
- Water and Air Syringe
- 750 ml Water Bottle
- Oil Compressor
- Vet Package
- Handle on Compressor

Technical Specifications	
Unit Dimensions	
Height	19 inches
Width	11 inches
Length	18 inches
Compressor	
Tank	2 gallons
Capacity	1/2 HP
Motor	110V
Flow Rates	4.1 CFM at 80 PSI
Weight	
Unit	65 pounds
Warranty	
	1 year



Item Number	Description
DTP00504	Countertop Unit with High-speed Push Button Handpiece, Slow-speed, Piezo
DTP00504-LED	Countertop Unit with High-speed LED Handpiece, Slow-speed, Piezo