

Countertop Dental Station



DTP00504



The DentalAire Countertop Unit is compact, powerful, and will easily integrate into your veterinary dentistry practice.

This countertop unit includes:

- High-speed Push Button Drill Handpiece
- Low-speed Polisher Handpiece
- LM Piezo Scaler with one #39 Tip Universal and #32 Flat Tip
- Water and Air Syringe
- 1 Liter Water Bottle
- Vet Package
- Handle on Compressor

Technical Specifications

Unit Dimensions

Height	19 inches
Width	11 inches
Length	18 inches

Compressor

Tank	2 gallons
Capacity	1/2 HP
Motor	110V
Flow Rates	4.1 CFM at 80 PSI

Weight

Unit	75 pounds
Shipping	100 pounds

Warranty

1 year

Item Number
DTP00504

Description
Countertop Unit with HS.SS.Piezo and Oil Compressor