

Diplomate Bur Kit

The DentalAire Bur Kit is an elegant solution for organizing your dental burs. The set includes 20 of the most common burs used in veterinary dentistry, all in an autoclavable bur block.

This bur kit unit includes:

- Round burs
- Cross-cut tapered fissure burs
- Conical white stones
- Two sizes of cylindrical diamond burs
- 12 extra sections for any additional burs



DTP10206

Included Burs

#1/2 Round	DTP10046	<ul style="list-style-type: none"> • Useful for making gutters (moats) around broken roots for retrieval. 	#700 Fissure	DTP10026	<ul style="list-style-type: none"> • For sectioning smaller and larger teeth
#1 Round	DTP10047	<ul style="list-style-type: none"> • For lateral bone removal during surgical extractions. • Use larger burs for larger teeth. 	#700L Fissure	DTP10030	<ul style="list-style-type: none"> • For sectioning smaller and larger teeth
#2 Round	DTP10048	<ul style="list-style-type: none"> • For lateral bone removal during surgical extractions. • Use larger burs for larger teeth. 	Feline Cylindrical	DTP835-014	<ul style="list-style-type: none"> • Useful for smoothing bone and retrieving broken root tips
#4 Round	DTP10052	<ul style="list-style-type: none"> • For lateral bone removal during surgical extractions. • Use larger burs for larger teeth. 	Canine Cylindrical	DTP835-018	<ul style="list-style-type: none"> • Useful for smoothing bone and retrieving broken root tips
#6 Round	DTP10054	<ul style="list-style-type: none"> • For lateral bone removal during surgical extractions. • Use larger burs for larger teeth. 	White Stone	DTP0450	<ul style="list-style-type: none"> • Ideal for conservatively smoothing bone • Used in the bonded sealant procedure