

Diplomate Wall Mount Dental Station



- LED high-speed handpiece, low-speed handpiece, LED piezo Ultrasonic Scaler
- Side mounted control adjustments
- 3-way syringe with quick change autoclavable tip
- Pressure adjustment for each handpiece
- Asepsis flex arm system with pneumatic brake, 52: reach, 27: vertical travel
- Oil Compressor - DA50-9D
- Soft touch variable speed foot control
- Dental Vet Package
- 7' umbilical tubing
- Wall mount plate included
- Dual 2 liter water bottles with wall mount or post mount option
- Nitrogen or compressor air driven
- Professional installation recommended (contact DentalAire)
- 2 - Year Warranty

Item Number

DTPJ-3300

DTPJ-3300-OL

DTPJ-3300-NC

Description

Diplomate Wall Mount Standard Style with Oil Compressor

Diplomate Wall Mount Unit with Oil-less Compressor

Diplomate Wall Mount Unit with No Compressor