

# DTX DIGITAL SENSORS

Instant results. Faster diagnosis.




Size 1 Sensor - DTP70008	
External Dimensions	25 x 39 mm
Active Image Area	20 x 30 mm
Pixel Size	20 x 20 μm
Number of Pixels	1.5 million
Theoretical Resolution	20 Lp/mm
Cable Length	2.5 meters



Size 2 Sensor - DTP70009	
External Dimensions	30.5 x 42 mm
External Dimensions with Boot	33 x 45 mm
Active Image Area	26 x 34 mm
Pixel Size	20 x 20 μm
Number of Pixels	2.2 million
Theoretical Resolution	20 Lp/mm
Cable Length	2.5 meters



Size 4c Sensor - DTP700020	
External Dimensions	42.1 x 57.9 mm
External Dimensions with Boot	N/A
Active Image Area	51.6 x 34 mm
Pixel Size	20 x 20 μm
Number of Pixels	2580 x 17.06
Cable Length	2.5 meters



■ Size 1   
 ■ Size 2   
 ■ Size 4c

## Contact Us

(800) 866-6881

[info@dentalaire.com](mailto:info@dentalaire.com)

3 year warranty

Sensor Replacement Program