

DTX Dental CR Reader

The Complete Image: CR Reader, DTX software, and imaging plates.

Maximize your practice's efficiency with the **DentalAire DTX CR Reader**, designed for rapid, high-definition processing of reusable imaging plates. Paired with our intuitive DTX software, you can seamlessly capture images, automate distribution, and generate comprehensive patient reports in seconds.

HIGH-QUALITY IMAGES

- Rapid, high-quality reading of reusable dental imaging plates
- Software optimizes acquisition, processing, and management
- Autocalibration ensures stability and consistency

FEATURES

- Automated, magnetized feeding tray
- "Push and Go" functionality
- Color LCD Touchscreen
- Intuitive user interface
- Instant image preview for quick quality verification
- Built-in erase function

DIMENSIONS

- Height: 10.4 in
- Width: 4.7 in
- Depth: 12.5 in
- Weight: 12.1 lbs



DTP710026



Unlimited Customer and Technical Support

TECHNICAL SPECIFICATIONS

CR Reader			
Selectable Pixel Size	High: 35µm Standard: 64 µm	Interface	USB 2.0 high-speed (480 Mbps) Ethernet (100 Mbps)
Sampling Pixel Pitch (35 µm)	14.3 lp/mm	Power Requirements	100 ~ 240V / 50 ~ 60Hz, 30W UPS required
Readout Time	4.1 ~ 7.2 seconds	Resolution	Data Capture: 16 bits per pixel, 65,000 graytones
Eraser Dimensions (H x W x D)	10.4 x 4.7 x 12.5		
Weight	12.1 pounds		
Warranty	1 year CR Reader		