


DTX Sensors



- Sensor 1**
- Sensor 2**
- Sensor 4c**



	Sensor 1	Sensor 2	Sensor 4c
External Dimensions	25 mm x 39 mm	30.5 mm x 42mm	42.1 mm x 57.9 mm
Active Image Area	20 mm x 30mm	26 mm x 34 mm	51.6mm x 34 mm
Pixel Size	20um x 20um	20um x 20um	20um x 20um
Number of Pixels	1.5 million	2.2 million	4.4 million
Cable Length	2.5 m	2.5 m	2.5 m
Resolution	20 lp/mm typ.	20 lp/mm typ.	20 lp/mm typ.
Warranty	3 years	3 years	3 years