



DENTALAIRE

Meeting and Exceeding the Needs of Veterinary Dentistry since 1988

DentalAire Diplomate Bur Kit

The DentalAire Diplomate Bur Set is an elegant solution for organizing your dental burs. The set includes 20 of the most common burs used in veterinary dentistry, all in an autoclavable bur block.



A selection of round burs, cross-cut tapered fissure burs, conical white stones and two sizes of cylindrical diamond burs are all you need for the majority of your dental procedures. Featuring 12 extra sections for your personal favorites.

Suggested uses for the burs are as follows:

- 1/2 Round- useful for making gutters (moats) around broken roots for retrieval.
- 1, 2, 4, 6 Round- for lateral bone removal during surgical extractions. Use larger burs for larger teeth.
- 700 and 701L Crosscut Tapered Fissure- for sectioning smaller and larger teeth respectively.
- Conical White Stones- ideal for conservatively smoothing bone and used in the bonded sealant procedure.
- Cylindrical Diamond Burs- useful for smoothing bone and retrieving broken root tips.



www.DentalAire.com | 800.866.6881

DentalAire | 17742 Mitchell North, Suite B Irvine CA 92614

info@dentalaire.com