

Dentalaire Diplomate Basic

Wall Mount Dental Station



Diplomate Basic

The **Dentaire Diplomate Basic** is a wall mount dental station that provides flexibility and efficiency for veterinary practices in need of space-saving solutions.

This newly designed dental station comes with three handpieces, including the *NEW* LED Swivel High-speed Handpiece. For easier access to internal components, this delivery systems includes the innovative flip top design.

- LED swivel high-speed handpiece
- Low-speed handpiece
- LED piezo ultrasonic scaler
- Air/water syringe
- Oil compressor
- Variable speed foot control
- One (1) - 2 liter water bottle
- 1-year warranty

DTP00705 Diplomate Basic
DTP00705-NC Diplomate Basic No Compressor



Customizable Scissor Arms

This dental station has mix and match arm options.

- Standard option includes two (2) 12 inch arms
- 6", 12", or 16" arm options
- Available in scissor configuration or single arm



High-speed LED Swivel Handpiece

The High-speed Swivel LED Handpiece provides direct light on the face of the handpiece.

- Eliminates fiber optic tubing
- Smooth and quiet operation



Provides flexibility and efficiency.


Included with the Diplomate Basic is the Dentalaire Vet Pack, which is a great starter kit for using your dental station. This kit includes:

- Mint prophy paste
- Snap-type prophy cups
- Snap-on prophy angle eco.
- Bur pack
- O-ring for tank drain
- All-In-One handpiece cleaner and lubricant



Looking to customize your dental station?

Contact us to learn more about your options.



Oil Compressor

- Horse Power: 1/2
- Noise Level: 40 decibels
- Tank Size: 2.25 gallons



High-quality Handpieces

- LED Swivel High-speed
- Low-speed
- LED Ultrasonic Scaler