

# Dentalaire Diplomate Standard

*Wall Mount Dental Station*



# Diplomate Standard

The **Dentalaire Diplomate Standard Style** is one of our most popular and versatile wall-mount dental systems. With side mounted control adjustments and a flex arm system with a pneumatic brake, this dental station can be designed to meet the exact specifications that you need.

This newly designed dental station comes with three handpieces, including the *NEW* LED Swivel High-speed Handpiece. For easier access to internal components, this delivery systems includes the innovative flip top design.

- LED swivel high-speed handpiece, low-speed handpiece, LED piezo ultrasonic scaler
- Air/water syringe
- Nitrogen or compressor air driven
- Variable speed foot control
- 2 - 2 liter water bottle
- 2 - year warranty excluding handpieces and accessories
- All other warranties are described in the individual product description



DTPJ-3300	Diplomate Wall Mount Standard Style with Oil Compressor
DTPJ-3300-OL	Diplomate Wall Mount Standard Style with Oil-less Compressor
DTPJ-3300-NC	Diplomate Wall Mount Standard Style with No Compressor
DTPJ-3300-P	Diplomate Wall Mount Standard Style with Plumbed Air/Water Regulators



## Asepsis Flex Arm System

The asepsis flex arm system with pneumatic brake provides more flexibility and efficiency .

- Reach: 52 inches
- Vertical Travel: 27 inches



## High-speed LED Swivel Handpiece

The High-speed Swivel LED Handpiece provides direct light on the face of the handpiece.

- Eliminates fiber optic tubing
- Smooth and quiet operation



## Dual 2-Liter Water Bottles


The water bottles can be placed in a highly convenient option that works best for your clinic.

- Wall or post mount option

## New wall-mount dental station.

Included with the Diplomat Standard is the Dentalaire Vet Pack, which is a great starter kit for using your dental station. This kit includes:

- Mint prophy paste
- Snap-type prophy cups
- Snap-on prophy angle eco.
- Bur pack
- O-ring for tank drain
- All-In-One handpiece cleaner and lubricant



### Oil Compressor

- Horse Power: 1/2
- Noise Level: 40 decibels
- Tank Size: 2.25 gallons



### High-quality Handpieces

- High-speed LED Swivel
- Low-speed
- Ultrasonic Scaler

**Looking to customize your dental station?**

Contact us to learn more about your options.