

Prestige Dental Station

VETERINARY MOBILE DENTAL STATION



Prestige Dental Station

Improve efficiency. Enhance patient care.

Enhance your clinic and improve dental procedures with **Dentalair's Prestige Dental Station**. As a top-choice mobile dental solution, this unit offers unmatched reliability, ample storage, and smooth operation for ultimate performance.

With the swivel LED high-speed handpiece, it is able to offer superior lighting that elevates your clinic's dental procedures.

HIGH-QUALITY HANDPIECES

- LED swivel high-speed
- Low-speed
- LED Piezo Scaler
- Air/water syringe

FEATURES

- Internal storage
- Asepsis, modern design with pneumatic lift
- Quick disconnect pressurized water connection
- Controls located on side panel
- Standard Midwest 4-hole connection
- 2 - 2 liter water bottle system
- Oil, oil-less, or no compressor option

DIMENSIONS

- Height: 39" - 30"
- Width: 22 1/2"
- Depth: 24" - 22"



✓ 2 year warranty *



OPTIMIZED FOR EFFICIENCY.

With your clinic's workflow in mind, the **Prestige Dental Station** gives you the freedom to easily store small tools and equipment into the large storage tray located under the lid.

In doing so, this can help provide your clinic the ample space it needs. Having an organized workspace can help maximize your productivity.

Lastly, repairs are easier and more efficient with parts accessible under the storage tray. Simply, open the lid, lift the storage tray, and easily access components for repair.


✓ Internal storage tray

✓ Enhanced serviceability

TAKE A CLOSER LOOK

Designed to help improve dental procedures, patient care, and clinic workflows, the Prestige Dental Station is equipped with innovative features to improve your workflow.

Explore how this mobile dental station can seamlessly integrate into your practice.



Asepsis Design

- Easily clean, disinfect, and create sterile environment
- User-friendly controls
- Smooth handpiece tubing



High-speed LED Swivel Handpiece

- Triple water spray cools entire surface of bur
- 410,000 RPM
- 360 °swivel rotation

TECHNICAL SPECIFICATIONS

Dentaire Item Number	Name
DTP00511	Prestige Dental Station with Oil Compressor
DTP00511-OL	Prestige Dental Station with Oil-less Compressor
DTP00512	Prestige Dental Station with No Compressor

Dental Station	
Height	39" - 30"
Width	22 1/2"
Depth	24" - 22"

Oil-less Compressor	
Horse Power	1
Noise Level	58 decibels
Tank Size	0.93 gallons

Oil Compressor	
Horse Power	1/2
Noise Level	40 decibels
Tank Size	2.25 gallons

Warranty	
Dental Station	2 years
Handpieces and Accessories	Warranties are described in individual product description



www.dentaire.com | 800.866.6881 | 17742 Mitchell North, STE B, Irvine, CA, 92614 | info@dentaire.com