

Prestige vs. Ultimate



	Prestige	Ultimate
Design	Asepsis, Modern	Traditional
Warranty	2 years	1 year
Piezo Scaler	LED Lighted Piezo Scaler	Non Lighted Piezo Scaler
Water Bottle System	2 – 2 Liter Water Bottle System	2 – 2 Liter Water Bottle System
Storage	Large Internal Storage Under Lid	No Internal Storage
Operation Controls	Side Located Recessed	Located on Front Panel
Service	Enhanced serviceability. All service parts are easily accessible under lid.	Delivery system with service parts located under top countertop
Height Adjustment	Pneumatic Lift	Manual
Compressor	Oil and Oil-less Available	Oil and Oil-less Available
High-speed Handpiece	LED Swivel	LED Swivel