

# Prestige or Ultimate Dental Station?

SOLUTIONS FOR EVERY CLINIC.

The **Prestige and Ultimate** mobile dental stations help improve veterinary dental procedures for clinics and ensure that each patient receives the highest quality of care they deserve.



Unlimited Customer and Technical Support



Included DentaLaire Dental Vet Package

## PRESTIGE DENTAL STATION



## ULTIMATE DENTAL STATION



<b>Design</b>	<b>Asepsis, Modern with Pneumatic Lift</b> <ul style="list-style-type: none"> <li>Helps promote a clean working station and eliminate harmful bacteria</li> </ul>	<b>Traditional with Manual Height Adjustment</b> <ul style="list-style-type: none"> <li>Ergonomic</li> </ul>
<b>Handpieces</b>	<b>LED Piezo Scaler</b> <ul style="list-style-type: none"> <li>LED Piezo Scaler, LED Swivel High-speed, Low-speed, Air/Water Syringe</li> </ul>	<b>Piezo Scaler, with No Light</b> <ul style="list-style-type: none"> <li>Piezo Scaler (with no light), LED Swivel High-speed, Low-speed, Air/Water Syringe</li> </ul>
<b>Water Bottle System</b>	<b>2 - 2 Liter Water Bottles</b> <ul style="list-style-type: none"> <li>Reliable water source during dental procedures and cleanings</li> </ul>	<b>2 - 2 Liter Water Bottles</b> <ul style="list-style-type: none"> <li>Reliable water source during dental procedures and cleanings</li> </ul>
<b>Operation Controls</b>	<b>Side Located Recessed</b> <ul style="list-style-type: none"> <li>Enhanced ergonomics, improved accessibility, and maximized workspace</li> </ul>	<b>Located on Front Panel</b> <ul style="list-style-type: none"> <li>Easy to access controls and user-friendly</li> </ul>
<b>Storage</b>	<b>Internal</b> <ul style="list-style-type: none"> <li>Large storage tray located under lid</li> </ul>	<b>No Storage</b> <ul style="list-style-type: none"> <li>Large workspace area that includes removable stainless steel tray</li> </ul>
<b>Service</b>	<b>Enhanced Serviceability</b> <ul style="list-style-type: none"> <li>All service parts are easily accessible under storage tray in lid</li> </ul>	<b>Service Parts Under Top Countertop</b> <ul style="list-style-type: none"> <li>Access service parts under top countertop</li> </ul>
<b>Compressor</b>	<b>Pick your option</b> <ul style="list-style-type: none"> <li>Oil, oil-less, and no compressor option</li> </ul>	<b>Pick your option</b> <ul style="list-style-type: none"> <li>Oil, oil-less, and no compressor option</li> </ul>
<b>Warranty</b>	<b>2 year warranty*</b> <ul style="list-style-type: none"> <li>Excludes handpieces and accessories</li> </ul>	<b>1 year warranty*</b> <ul style="list-style-type: none"> <li>Excludes handpieces and accessories</li> </ul>

