



DENTALAIRE

Meeting and Exceeding the Needs of Veterinary Dentistry since 1988

Ram Micro Motor

MICROTORQUE DIGITAL CONTROL BOX FEATURES:

- 0-35,000 RPM
- Digital Display for Precise Speed Control
- Foot Control
- Standby Mode for Handpiece Motor
- Forward or Reverse Rotation
- Microprocessor Controlled for Safety and Reliability
- Microprocessor for Audible Signals, Soft Start and Stop
- Microprocessor with Memory Functions
- Improved Overload Sensor Feedback
- Touch Panel Buttons for Durability
- No Ventilation Needed to Prolong Life by Keeping Out Dust and Moisture
- Improved Speed Control and Variation

SLIM NOSE HANDPIECE FEATURES:

- 35,000 RPM
- 3/32" Chuck
- Tapered Slim Style Nose Tip
- Twist Type Chuck Release



DTP103RMK

MICROTORQUE DIGITAL CONTROL BOX SPECIFICATIONS:

- Input: AC110V/230V 50HZ-60HZ
- Output: DC 0-32V
- Max Amp: 2A
- Rotate Speed: 0-35,000rpm
- Power: 60W
- Dimension: 6" Length x 3 1/4" Width x 5 1/4" Height
- Weight: 4.4lbs.
- 6.75" Length x 4" Width x 2.75" Height

VARIABLE FOOT PEDAL SPECIFICATIONS:

- 6.75" Length x 4" Width x 2.75" Height

SLIM NOSE HANDPIECE SPECIFICATIONS:

- Torque: 3.2Ncm
- Size: 6" x 1.12"
- Net Weight: 8oz.

INCLUDED VET PACKAGE:

- Snap-on Prophy Angle Prophy Paste (Fruit Flavored)
- Snap-type Prophy Cups (Pack of 144)
- Latch Type Contra Angle Eco.
- Dia. Disc Wheel Double Side/Mand
- Latch Burs As. (Pack of 5)

www.DentalAireProducts.com | 800.866.6881

DentalAire Products | 17742 Mitchell North, Suite B Irvine CA 92614

info@dentalAireproducts.com