

Dentalaire Ultimate Dental Station

Mobile Dental Station



Ultimate Dental Station

The **Ultimate Dental Station** has a large workspace area that includes a removable stainless-steel tray, which can be used to place dental instruments and supplies.

The included handpieces offer versatility with the new LED Swivel high-speed handpiece, built-in piezo scaler, and a low-speed handpiece.

- Telescopic height adjustment from 40" down to 26.5"
- LED Swivel high-speed handpiece
- Low-speed handpiece
- Piezo scaler
- Three-way syringe
- Stainless steel tray
- Quick disconnect pressurized water connection
- 1- year warranty excluding handpieces and accessories
- All other warranties are described in the individual product description

DTP00509 Ultimate Dental Station with Oil Compressor

DTP00509-OL Ultimate Dental Station with Oil-less Compressor



NEW High-speed LED Swivel Handpiece

The High-speed Swivel LED Handpiece's macro-angled design and high color rendering index provides direct light on the face of the handpiece that is high in intensity and provides larger coverage.

- Eliminates fiber optic tubing
- Swivel rotation provides smooth and quiet operation
- Helps reduce arm fatigue



Large Work Area

The stainless steel removable tray can be used to place dental instruments and supplies.

- 22" x 19" work area
- Removable stainless steel tray

Complete with everything you need.

Included with the Ultimate Dental Stations is the Dentalaire Vet Pack, which is a great starter kit for using your dental station. This kit includes:

- Mint prophy paste
- Snap-type prophy cups
- Snap-on prophy angle eco.
- Bur pack
- O-ring for tank drain
- All-In-One handpiece cleaner and lubricant



Looking to customize your dental station?
Contact us to learn more about your options.



Oil Compressor

- Horse Power: 1/2
- Noise Level: 40 decibels
- Tank Size: 2.25 gallons



Adjustable Height

- 40" down to 26.5"
- Dimensions: 40 -26.5" H x 22" D x 22" W



High-quality Handpieces

- High-speed LED Swivel
- Low-speed
- Ultrasonic Scaler